

As an Interpath customer who receives electronic results or sends electronic orders you may need to be notified when we update our Service Manual. Although we try to keep these changes to a minimum, laboratory medicine is an evolving industry requiring changes to our technology from time to time. Depending on the requirements of your EMR or Hospital Information System you may be required to make similar changes to your system in order to correctly process inbound electronic results and create outbound electronic orders.

If you are uncertain that you are required to update your system we recommend that you contact your vendor for more information. As your laboratory service provider we are available to participate in the discussion with your vendor so that you clearly understand the impact of these changes.

Included in this email:

- This cover letter with a summary of the changes
- Microsoft Word® Document with the detail of these changes to our Service Manual
- Interpath Master Order/Result Compendium

Additional information including our most recent Service Manual and additional contact information can be found at www.interpathlab.com

Effective Date: February 16th, 2021

Order Code	Test Name	NC	CC	CPT	SRC	RRC	NT	DT	AOE
		Name Change	Component Change	CPT Change	Specimen Requirements Change	Reference Range Change	New Test	Discontinued Test	Ask on Order Entry Questions
93627	Fluphenazine				◆	◆			
91345	Haloperidol					◆			
90200	Hemoglobin Evaluation Reflexive Cascade		◆						
93092	IgG Subclass 4				◆	◆			
91143	Immunoglobulin G Subclasses				◆	◆			
91609	Interleukin 28 B (IL28B)- Associated Variants, 2 SNPs							◆	
90300	Lead, Capillary							◆	
76043	Lead, Capillary						◆		
91281	Maternal AFP		◆						◆
91011	Maternal Quad Panel		◆						◆
93638	Maternal Screen #1		◆						◆
90070	Meprobamate							◆	
91174	Oligoclonal Band Profile				◆	◆			
91293	Phosphatidylserine Abs, IgGAM								
2739	Procalcitonin						◆		
93845	Procalcitonin							◆	

93627 Fluphenazine

SRC/RRC

Specimen:	
Collect:	One Red Top Also Acceptable One Lavender (EDTA) One Pink Top (EDTA)
Submit:	1 mL (Min:0.5 mL) Serum. Submit Refrigerated. Submit in a Standard Transport Tube. Also Acceptable 1 mL (Min:1 mL) Plasma. Submit Refrigerated. Submit in a Standard Transport Tube.
Special Handling:	Separate from cells ASAP Timing of specimen collection: Pre-dose (trough) draw - At steady state concentration.
Rejection Criteria:	Hemolyzed specimens Use of separator tubes Whole blood
Stability:	Ambient: 2 Day(s); Refrigerated: 1 Week(s); Frozen: 1 Month(s); Incubated: Unacceptable
Methodology:	Quantitative Liquid Chromatography-Tandem Mass Spectrometry
Performed:	Monday, Wednesday, Friday
Reported:	2-6 Day(s)
CPT Codes:	80342

Please take note of changes to rejection criteria and reference range.

Reference Range:

Therapeutic Range: 1.0-10.0 ng/mL

Toxic: >15 ng/mL

91345 Haloperidol

RRC

Please take note of change to reference range.

Reference Range:

Therapeutic Range: 5.0-20.0 ng/mL

Toxic: >50 ng/mL

90200 Hemoglobin Evaluation Reflexive Cascade
CC

Specimen:																			
Collect:	One Lavender (EDTA) Also Acceptable One Pink Top (EDTA)																		
Submit:	4 mL (Min:2 mL) Whole blood in Lavender (EDTA). Submit Refrigerated. Also Acceptable 4 mL (Min:2 mL) Whole blood in Pink Top (EDTA). Submit Refrigerated.																		
Special Handling:	Patient history form, including information from a recent CBC, is required for interpretation. Separate specimens must be submitted when multiple tests are ordered. Abnormal results reflex. Additional charges will apply.																		
Stability:	Ambient: Unacceptable; Refrigerated: 1 Week(s); Frozen: Unacceptable; Incubated: Unacceptable																		
Methodology:	Electrophoresis; Fluorescence Resonance Energy Transfer; High Performance Liquid Chromatography; Polymerase Chain Reaction (PCR); RBC Solubility; Sequencing																		
Performed:	Sun-Sat																		
Reported:	Varies																		
CPT Codes:	83021																		
Interpretive Data:	Please see report for interpretive data.																		
Components:	<table border="0"> <tr> <td>93437 - HEMOGLOBIN A2</td> <td>93438 - HEMOGLOBIN F</td> </tr> <tr> <td>93439 - HEMOGLOBIN S</td> <td>93440 - HEMOGLOBIN C</td> </tr> <tr> <td>93441 - HEMOGLOBIN E</td> <td>93442 - HEMOGLOBIN OTHER</td> </tr> <tr> <td>93436 - HEMOGLOBIN A1</td> <td>93443 - HEMOGLOBIN INTERP</td> </tr> <tr> <td>93619 - SICKLE CELL SOL</td> <td>93977 - HGB, CAP ELECRO</td> </tr> <tr> <td>90187 - ALPHA GLOBIN</td> <td>90188 - ALPHA THALASSEMIA</td> </tr> <tr> <td>90199 - BETA GLOBIN SEQU</td> <td>92019 - BETA GLOBIN DEL</td> </tr> <tr> <td>92026 - HEMOGLOBIN LEPORE</td> <td>92027 - INTERPRETATION</td> </tr> <tr> <td>92182 - GAMMA GLOBIN SEQ</td> <td></td> </tr> </table>	93437 - HEMOGLOBIN A2	93438 - HEMOGLOBIN F	93439 - HEMOGLOBIN S	93440 - HEMOGLOBIN C	93441 - HEMOGLOBIN E	93442 - HEMOGLOBIN OTHER	93436 - HEMOGLOBIN A1	93443 - HEMOGLOBIN INTERP	93619 - SICKLE CELL SOL	93977 - HGB, CAP ELECRO	90187 - ALPHA GLOBIN	90188 - ALPHA THALASSEMIA	90199 - BETA GLOBIN SEQU	92019 - BETA GLOBIN DEL	92026 - HEMOGLOBIN LEPORE	92027 - INTERPRETATION	92182 - GAMMA GLOBIN SEQ	
93437 - HEMOGLOBIN A2	93438 - HEMOGLOBIN F																		
93439 - HEMOGLOBIN S	93440 - HEMOGLOBIN C																		
93441 - HEMOGLOBIN E	93442 - HEMOGLOBIN OTHER																		
93436 - HEMOGLOBIN A1	93443 - HEMOGLOBIN INTERP																		
93619 - SICKLE CELL SOL	93977 - HGB, CAP ELECRO																		
90187 - ALPHA GLOBIN	90188 - ALPHA THALASSEMIA																		
90199 - BETA GLOBIN SEQU	92019 - BETA GLOBIN DEL																		
92026 - HEMOGLOBIN LEPORE	92027 - INTERPRETATION																		
92182 - GAMMA GLOBIN SEQ																			

Please take note of change to components.
Component Changes:
Remove: 92026 - HEMOGLOBIN LEPORE
Add: 92182 - GAMMA GLOBIN SEQ

93092 IgG Subclass 4
SRC/RRC

Specimen:	
Collect:	One SST Also Acceptable One Red Top
Submit:	1 mL (Min:0.5 mL) Serum. Submit Refrigerated. Submit in a Standard Transport Tube.
Rejection Criteria:	Grossly Hemolyzed Samples Lipemic Samples
Stability:	Ambient: Unacceptable; Refrigerated: 2 Week(s); Frozen: 6 Month(s); Incubated: Unacceptable
Methodology:	Quantitative Immunoturbidimetry
Performed:	Sun-Sat
Reported:	2-4 Day(s)
CPT Codes:	82787

Please take note of changes to rejection criteria, stability, methodology, performed dates, and reference ranges.

Reference Ranges:

Age	Reference Range
0-2 years	1-34 mg/dL
3-4 years	1-65 mg/dL
5-9 years	0-168 mg/dL
10-14 years	1-103 mg/dL
15-18 years	2-170 mg/dL
19 and older	1-123 mg/dL

91143 Immunoglobulin G Subclasses
SRC/RRC

Specimen:	
Collect:	One SST Also Acceptable One Red Top
Submit:	2 mL (Min:1 mL) Serum. Submit Refrigerated. Submit in a Standard Transport Tube.
Special Handling:	Separate from cells ASAP
Rejection Criteria:	Grossly Hemolyzed Samples Lipemic Samples
Stability:	Ambient: Unacceptable; Refrigerated: 2 Week(s); Frozen: 6 Month(s); Incubated: Unacceptable
Methodology:	Quantitative Immunoturbidimetry
Performed:	Sun-Sat
Reported:	2-4 Day(s)
CPT Codes:	82787x4
Interpretive Data:	Please see report for interpretive data.
Components:	93089 - IgG SUBCLASS 1 93090 - IgG SUBCLASS 2 93091 - IgG SUBCLASS 3 93092 - IgG SUBCLASS 4

Please take note of changes to rejection criteria, stability, methodology, performed dates, and reference ranges.

Reference Ranges:
93089 – IgG Subclass 1

Age	Reference Range
0-2 years	167-900 mg/dL
3-4 years	313-941 mg/dL
5-9 years	363-1276 mg/dL
10-14 years	316-1076 mg/dL
15-18 years	325-894 mg/dL
19 and older	240-1118 mg/dL

93090 – IgG Subclass 2

Age	Reference Range
0-2 years	55-359 mg/dL
3-4 years	72-287 mg/dL
5-9 years	27-398 mg/dL
10-14 years	86-509 mg/dL
15-18 years	156-625 mg/dL
19 and older	124-549 mg/dL

Changes continued on following page

93091 – IgG Subclass 3

Age	Reference Range
0-2 years	34-85 mg/dL
3-4 years	25-117 mg/dL
5-9 years	17-169 mg/dL
10-14 years	14-201 mg/dL
15-18 years	34-246 mg/dL
19 and older	21-134 mg/dL

93092 – IgG Subclass 4

Age	Reference Range
0-2 years	1-34 mg/dL
3-4 years	1-65 mg/dL
5-9 years	0-168 mg/dL
10-14 years	1-103 mg/dL
15-18 years	2-170 mg/dL
19 and older	1-123 mg/dL

91609 Interleukin 28 B (IL28B)- Associated Variants, 2 SNPs
Test is being discontinued.

DT

90300 Lead, Capillary
DT

Please take note test is being discontinued and replaced with Test 76043 Lead, Capillary.

76043 Lead, Capillary
NT

Specimen:	
Collect:	One Lavender (EDTA)
Submit:	0.5 mL (Min:0.3 mL) Whole blood in Lavender (EDTA). Submit Ambient.
Special Handling:	Invert 10 times to prevent specimen clotting. Stability: If the specimen is drawn and stored in the appropriate container the lead values do not change with time.
Rejection Criteria:	Specimens collected/transported in tubes other than Lavender Pediatric (EDTA) or Trace Element Free Transport Tube. Heparin anticoagulant or clotted specimens.
Stability:	Ambient: 24 Month(s); Refrigerated: 24 Month(s); Frozen: Unacceptable; Incubated: Unacceptable
Methodology:	Inductively Coupled Plasma Mass Spectrometry
Performed:	Mon-Fri
Reported:	2-8 Day(s)
CPT Codes:	83655
Interpretive Data:	General Reference Range : <=4.9 ug/dL

New test available for order.

91281 Maternal AFP
CC/AOE

Specimen:																	
Collect:	One SST Also Acceptable One Red Top																
Submit:	1 mL (Min:0.5 mL) Serum. Submit Refrigerated. Submit in a Standard Transport Tube.																
Special Handling:	Avoid Repeated Freeze/Thaw Cycles Specimen must be drawn between 14 weeks, 0 days and 24 weeks, 6 days gestation. Special Form - Patient History for Maternal Serum Testing The following information is required and must accompany the sample for test Interpretation: Patient's date of birth, current weight, due date, dating method (US, LMP), number of fetuses present, patient's race, if the patient was diabetic at the time of conception, if there is a known family history of neural tube defects, if the patient is currently smoking, if the patient is taking valproic acid or carbamazepine (Tegretol), and if this is a repeat sample.																
Rejection Criteria:	Grossly Hemolyzed Samples Plasma																
Stability:	Ambient: 3 Day(s); Refrigerated: 2 Week(s); Frozen: 12 Month(s); Incubated: Unacceptable																
Methodology:	Quantitative Chemiluminescent Immunoassay																
Performed:	Sun-Sat																
Reported:	3-4 Day(s)																
CPT Codes:	82105																
Interpretive Data:	Please see report for interpretive data.																
Components:	<table border="0"> <tr> <td>93241 - SPECIMEN</td> <td>93248 - Maternal Age</td> </tr> <tr> <td>93242 - DATING</td> <td>93202 - Est. Due Date</td> </tr> <tr> <td>93252 - Gestational Age</td> <td>93238 - Maternal Weight</td> </tr> <tr> <td>93243 - Mat. Insulin DM</td> <td>93244 - Family Hx of NTD</td> </tr> <tr> <td>93245 - Maternal Race</td> <td>93246 - Number of Fetuses</td> </tr> <tr> <td>93169 - Patients AFP</td> <td>93239 - MoM for AFP</td> </tr> <tr> <td>93240 - Interpretation</td> <td>92052 - SMOKING</td> </tr> <tr> <td>92184 - In Vitro Fert</td> <td>92185 - Donor egg age</td> </tr> </table>	93241 - SPECIMEN	93248 - Maternal Age	93242 - DATING	93202 - Est. Due Date	93252 - Gestational Age	93238 - Maternal Weight	93243 - Mat. Insulin DM	93244 - Family Hx of NTD	93245 - Maternal Race	93246 - Number of Fetuses	93169 - Patients AFP	93239 - MoM for AFP	93240 - Interpretation	92052 - SMOKING	92184 - In Vitro Fert	92185 - Donor egg age
93241 - SPECIMEN	93248 - Maternal Age																
93242 - DATING	93202 - Est. Due Date																
93252 - Gestational Age	93238 - Maternal Weight																
93243 - Mat. Insulin DM	93244 - Family Hx of NTD																
93245 - Maternal Race	93246 - Number of Fetuses																
93169 - Patients AFP	93239 - MoM for AFP																
93240 - Interpretation	92052 - SMOKING																
92184 - In Vitro Fert	92185 - Donor egg age																

Please take note of changes to components.
Component Changes:

Add: 92184 - In Vitro Fert; 92185 - Donor egg age

AOE changes:

Add: Q1147 – In Vitro Fertilization; Q1165 – Donor Egg Age

91011 Maternal Quad Panel

CC/AOE

Specimen:																										
Collect:	One SST Also Acceptable One Red Top																									
Submit:	3 mL (Min:1 mL) Serum. Submit Refrigerated. Submit in a Standard Transport Tube.																									
Special Handling:	Avoid Repeated Freeze/Thaw Cycles Separate aliquot required for each frozen test ordered Separate from cells ASAP Collect at 14-24 weeks gestation only. Remarks: The following information is required and must accompany the sample for test interpretation: Patient's date of birth, current weight, due date, dating method (US, LMP), number of fetuses present, patient's race, if the patient was diabetic at the time of conception, if there is a known family history of neural tube defects, if the patient has had a previous pregnancy with a trisomy, if the patient is currently smoking, if the patient is taking valproic acid or carbamazepine (Tegretol), if this is a repeat sample, and the age of the egg donor if in vitro fertilization.																									
Rejection Criteria:	Grossly Hemolyzed Samples Plasma																									
Stability:	Ambient: 3 Day(s); Refrigerated: 2 Week(s); Frozen: 12 Month(s); Incubated: Unacceptable																									
Methodology:	Quantitative Chemiluminescent Immunoassay																									
Performed:	Sun-Sat																									
Reported:	3-5 Day(s)																									
CPT Codes:	82105 84702	82677 86336																								
Interpretive Data:	Please see report for interpretive data.																									
Components:	<table border="0"> <tr> <td>93169 - Patients AFP</td> <td>93202 - Est. Due Date</td> </tr> <tr> <td>93206 - MoM FOR DIA</td> <td>93207 - Patients DIA</td> </tr> <tr> <td>93238 - Maternal Weight</td> <td>93239 - MoM for AFP</td> </tr> <tr> <td>93240 - Interpretation</td> <td>93241 - SPECIMEN</td> </tr> <tr> <td>93242 - DATING</td> <td>93243 - Mat. Insulin DM</td> </tr> <tr> <td>93244 - Family Hx of NTD</td> <td>93245 - Maternal Race</td> </tr> <tr> <td>93246 - Number of Fetuses</td> <td>93248 - Maternal Age</td> </tr> <tr> <td>93249 - Patients hCG</td> <td>93251 - MoM for hCG</td> </tr> <tr> <td>93252 - Gestational Age</td> <td>93253 - Patients uE3</td> </tr> <tr> <td>93255 - MoM for uE3</td> <td>91593 - Hx of Aneuploidy</td> </tr> <tr> <td>90059 - EER SCREEN</td> <td>92052 - SMOKING</td> </tr> <tr> <td>92185 - Donor egg age</td> <td></td> </tr> </table>		93169 - Patients AFP	93202 - Est. Due Date	93206 - MoM FOR DIA	93207 - Patients DIA	93238 - Maternal Weight	93239 - MoM for AFP	93240 - Interpretation	93241 - SPECIMEN	93242 - DATING	93243 - Mat. Insulin DM	93244 - Family Hx of NTD	93245 - Maternal Race	93246 - Number of Fetuses	93248 - Maternal Age	93249 - Patients hCG	93251 - MoM for hCG	93252 - Gestational Age	93253 - Patients uE3	93255 - MoM for uE3	91593 - Hx of Aneuploidy	90059 - EER SCREEN	92052 - SMOKING	92185 - Donor egg age	
93169 - Patients AFP	93202 - Est. Due Date																									
93206 - MoM FOR DIA	93207 - Patients DIA																									
93238 - Maternal Weight	93239 - MoM for AFP																									
93240 - Interpretation	93241 - SPECIMEN																									
93242 - DATING	93243 - Mat. Insulin DM																									
93244 - Family Hx of NTD	93245 - Maternal Race																									
93246 - Number of Fetuses	93248 - Maternal Age																									
93249 - Patients hCG	93251 - MoM for hCG																									
93252 - Gestational Age	93253 - Patients uE3																									
93255 - MoM for uE3	91593 - Hx of Aneuploidy																									
90059 - EER SCREEN	92052 - SMOKING																									
92185 - Donor egg age																										

Please take note of changes to components.

Component Changes:

Add: 92185 - Donor egg age

AOE changes:

Add: Q1165 – Donor Egg Age

93638 Maternal Screen #1
CC/AOE

Specimen:																							
Collect:	One SST Also Acceptable One Red Top																						
Submit:	0.5 mL (Min:0.3 mL) Serum. Submit Refrigerated. Submit in a Standard Transport Tube.																						
Special Handling:	Specimen must be drawn between 10 weeks, 0 days and 13 weeks, 6 days gestation (Crown-Rump length (CRL) must be 32.4-83.9 mm). Requires complete gestational information including CRL. This test also requires the following information: Patient's date of birth, current weight, number of fetuses present, patient's race, if the patient was diabetic at the time of conception, if there is a known family history of neural tube defects, if the patient has had a previous pregnancy with a trisomy, if the patient is currently smoking, if the patient is taking valproic acid or carbamazepine (Tegretol), if this is a repeat sample, and the age of the egg donor if in vitro fertilization.																						
Rejection Criteria:	Hemolyzed specimens Plasma Specimens exposed to repeat freeze/thaw cycles																						
Stability:	Ambient: 3 Day(s); Refrigerated: 2 Week(s); Frozen: 3 Month(s); Incubated: Unacceptable																						
Methodology:	Quantitative Chemiluminescent Immunoassay																						
Performed:	Sun-Sat																						
Reported:	3-5 Day(s) Final interpretative report available when second specimen testing is complete.																						
CPT Codes:	84163																						
Interpretive Data:	Please see report for interpretive data.																						
Components:	<table border="0"> <tr> <td>93963 - Patients PAPP-A</td> <td>93961 - Nuchal Trans (NT)</td> </tr> <tr> <td>93240 - Interpretation</td> <td>93248 - Maternal Age</td> </tr> <tr> <td>93238 - Maternal Weight</td> <td>93202 - Est. Due Date</td> </tr> <tr> <td>93252 - Gestational Age</td> <td>93246 - Number of Fetuses</td> </tr> <tr> <td>93245 - Maternal Race</td> <td>91593 - Hx of Aneuploidy</td> </tr> <tr> <td>93962 - Crown Rump Length</td> <td>93965 - Sonographer Cert</td> </tr> <tr> <td>93966 - Sonographer Name</td> <td>93967 - Ultrasound Date</td> </tr> <tr> <td>93969 - Date for Sample 2</td> <td>93970 - EER Maternal Scrn</td> </tr> <tr> <td>92053 - NUCHAL TWIN B</td> <td>93242 - DATING</td> </tr> <tr> <td>92052 - SMOKING</td> <td>93241 - SPECIMEN</td> </tr> <tr> <td>92054 - Crown Rump Twin B</td> <td>92185 - Donor egg age</td> </tr> </table>	93963 - Patients PAPP-A	93961 - Nuchal Trans (NT)	93240 - Interpretation	93248 - Maternal Age	93238 - Maternal Weight	93202 - Est. Due Date	93252 - Gestational Age	93246 - Number of Fetuses	93245 - Maternal Race	91593 - Hx of Aneuploidy	93962 - Crown Rump Length	93965 - Sonographer Cert	93966 - Sonographer Name	93967 - Ultrasound Date	93969 - Date for Sample 2	93970 - EER Maternal Scrn	92053 - NUCHAL TWIN B	93242 - DATING	92052 - SMOKING	93241 - SPECIMEN	92054 - Crown Rump Twin B	92185 - Donor egg age
93963 - Patients PAPP-A	93961 - Nuchal Trans (NT)																						
93240 - Interpretation	93248 - Maternal Age																						
93238 - Maternal Weight	93202 - Est. Due Date																						
93252 - Gestational Age	93246 - Number of Fetuses																						
93245 - Maternal Race	91593 - Hx of Aneuploidy																						
93962 - Crown Rump Length	93965 - Sonographer Cert																						
93966 - Sonographer Name	93967 - Ultrasound Date																						
93969 - Date for Sample 2	93970 - EER Maternal Scrn																						
92053 - NUCHAL TWIN B	93242 - DATING																						
92052 - SMOKING	93241 - SPECIMEN																						
92054 - Crown Rump Twin B	92185 - Donor egg age																						

Please take note of changes to componnets.
Component Changes:

Add: 92185 - Donor egg age

AOE changes:

Add: Q1165 – Donor Egg Age

90070 Meprobamate
Test is being discontinued.

DT

91174 Oligoclonal Band Profile

SRC/RRC

Specimen:	
Collect:	One SST CSF in Sterile Specimen Container Also Acceptable One Red Top
Submit:	1 mL (Min:0.6 mL) Serum. Submit Refrigerated. Submit in a Standard Transport Tube. 1.5 mL (Min:1 mL) CSF in Sterile Specimen Container. Submit Refrigerated.
Special Handling:	Allow specimen to clot completely at room temperature Separate from cells ASAP Serum specimen should be drawn within 48 hours of CSF collection
Rejection Criteria:	Grossly Hemolyzed Samples Grossly Lipemic Samples
Stability:	Ambient: Unacceptable; Refrigerated: 2 Week(s); Frozen: 6 Month(s); Incubated: Unacceptable
Methodology:	Electrophoresis; Qualitative Isoelectric Focusing; Quantitative Immunoturbidimetry
Performed:	Monday, Wednesday, Friday
Reported:	2-5 Day(s)
CPT Codes:	82040 82042 83873 83916 82784x2
Interpretive Data:	Please see report for interpretive data.
Components:	93082 - IgG, SERUM 93094 - IgG, CSF 93071 - ALBUMIN, CSF 93095 - ALBUMIN, SERUM 93096 - ALBUMIN INDEX 93097 - IgG SYNTHESIS 93191 - IgG INDEX 93190 - CSF IgG/ALB RATIO 93181 - OLIGO BANDS 93867 - OLIGO BANDS NUMBER 93182 - INTERPRETATION

Please take note of changes to submit volumes, rejection criteria, stability, methodology and reference range.

Reference Range:

93082 – IgG, SERUM

Age	Reference Range
0-2 years	242-1108 mg/dL
3-4 years	485-1160 mg/dL
5-9 years	514-1672 mg/dL
10-14 years	581-1652 mg/dL
15-18 years	479-1433 mg/dL
19 and older	768-1632 mg/dL

91293 Phosphatidylserine Abs, IgGAM
SRC/RRC

Specimen:	
Collect:	One SST Also Acceptable One Red Top
Submit:	0.5 mL (Min:0.3 mL) Serum. Submit Refrigerated. Submit in a Standard Transport Tube.
Special Handling:	Avoid Repeated Freeze/Thaw Cycles
Rejection Criteria:	Bacterially Contaminated Samples Heat inactivated Hemolyzed specimens Lipemic Samples
Stability:	Ambient: 2 Day(s); Refrigerated: 2 Week(s); Frozen: 1 Month(s) ; Incubated: Unacceptable
Methodology:	Semi-Quantitative Enzyme-Linked Immunosorbent Assay
Performed:	Tuesday, Wednesday, Friday, Saturday, Sunday
Reported:	2-4 Day(s)
CPT Codes:	86148
Interpretive Data:	Please see report for interpretive data.
Components:	93585 - PHOSPHATIDYL IgA 93587 - PHOSPHATIDYL IgG 93586 - PHOSPHATIDYL IgM

Please take note of changes to stability and reference range.

Reference Ranges:

93585- PHOSPHATIDYL IgA: <20 APS (IgA antiphosphatidylserine units)
 93586- PHOSPHATIDYL IgM: <22 MPS (IgM antiphosphatidylserine units)
 93587- PHOSPHATIDYL IgG: <16 GPS (IgG antiphosphatidylserine units)

2739 Procalcitonin
NT

Specimen:	
Collect:	One SST Also Acceptable One PPT White Top (EDTA)
Submit:	1 mL (Min:0.5 mL) Serum. Submit Refrigerated. Submit in a Standard Transport Tube. Also Acceptable 1 mL (Min:0.5 mL) Plasma. Submit Refrigerated. Submit in a Standard Transport Tube.
Special Handling:	Separate from cells within 2 hours of collection
Stability:	Ambient: 1 Day(s); Refrigerated: 2 Day(s); Frozen: 3 Month(s); Incubated: Unacceptable
Methodology:	Electrochemiluminescence Immunoassay (ECLIA)
Performed:	Sun-Sat
Reported:	1-2 Day(s)
CPT Codes:	84145
Interpretive Data:	<p>General Reference Range : 0.02-0.08 ng/mL</p> <p>Procalcitonin Interpretation: <0.5 ng/mL - Associated with low risk for progression to severe sepsis on the first day of ICU admission. This does not exclude an infection, especially with localized infections without systemic signs that could be associated with such low concentrations or in cases of initial stages of infection (<6hrs). 2.0 ng/mL - Associated with a high risk for progression to severe sepsis and/or septic shock.</p> <p>Biotin in specimens taken from patients on high-dose biotin therapy or supplements may interfere with this test and cause inaccurate test results. It is recommended that for patients receiving therapy with high biotin doses (> 5 mg/day), no laboratory test specimen should be collected until at least 8 hours after the last biotin administration.</p>

New test available for order.
93845 Procalcitonin
DC
Please take note test is being discontinued and replaced with Test# 2739 Procalcitonin.