



Interpath Laboratory, Inc. Test File Update

As an Interpath customer who receives electronic results or sends electronic orders you may need to be notified when we update our Service Manual. Although we try to keep these changes to a minimum, laboratory medicine is an evolving industry requiring changes to our technology from time to time. Depending on the requirements of your EMR or Hospital Information System you may be required to make similar changes to your system in order to correctly process inbound electronic results and create outbound electronic orders.

If you are uncertain that you are required to update your system we recommend that you contact your vendor for more information. As your laboratory service provider we are available to participate in the discussion with your vendor so that you clearly understand the impact of these changes.

Included in this email:

- This cover letter with a summary of changes
- Microsoft Word® Document with the detail of these changes to our Service Manual
- Interpath Master Order/Result Compendium

Additional information including our most recent Service Manual and additional contact information can be found at www.interpathlab.com

Effective Date: May 3, 2021

Order Code	Test Name	NC	CC	CPT	SRC	RRC	NT	DT	AOE
		Name Change	Component Change	CPT Change	Specimen Requirements Change	Reference Range Change	New Test	Discontinued Test	Ask on Order Entry Questions
70077	Amphetamine, Urine (Confirmation)							◆	
74020	Amphetamine, Urine Confirmation						◆		
70128	Antidepressant (SSRI), Urine (Confirmation)							◆	
74030	Antidepressant (SSRI), Urine Confirmation						◆		
74040	Antidepressant (SSRI 2), Urine Confirmation						◆		
74050	Antidepressant (TCA), Urine Confirmation						◆		
70080	Barbiturates, Urine (Confirmation)							◆	
74060	Barbiturates, Urine Confirmation						◆		
70089	Benzodiazepines, Urine (Confirmation)							◆	
74070	Benzodiazepine, Urine Confirmation						◆		
70086	Buprenorphine, Urine (Confirmation)							◆	
74080	Buprenorphine, Urine Confirmation						◆		
70133	Cannabinoids, Urine (Confirmation)							◆	
74017	Cannabinoids, Urine Confirmation						◆		
2759	Cardiovascular Disease Risk Panel + NMR, Lipoprofile		◆						
70101	Carisoprodol, Urine (Confirmation)							◆	
74090	Carisoprodol, Urine Confirmation						◆		
70100	Cocaine, Urine (Confirmation)							◆	
74062	Cocaine, Urine Confirmation						◆		
74100	Codeine + Morphine, Urine Confirmation						◆		
70104	Cyclobenzaprine, Urine (Confirmation)							◆	
74022	Cyclobenzaprine, Urine Confirmation						◆		
74110	Dextromethorphan, Urine Confirmation						◆		
74120	Doxepin, Urine Confirmation						◆		
6955	Drug Screen, Urine Medical							◆	
73010	Drug Screen, Urine Medication Management						◆		◆
73020	Drug Screen, Urine Medication Management w/o THC and Alcohol						◆		◆
7000	Drug Screen, Urine Pain Management							◆	
73000	Drug Screen, Urine Pain Management						◆		
70137	Ecstasy, Urine (Confirmation)							◆	
74130	Ecstasy, Urine Confirmation						◆		
74140	Embeda, Urine Confirmation						◆		
70105	Ethyl Alcohol, Urine (Confirmation)							◆	
74150	Ethyl Alcohol, Urine Confirmation						◆		

		NC	CC	CPT	SRC	RRC	NT	DT	AOE
70108	Fentanyl, Urine (Confirmation)							◆	
74160	Fentanyl, Urine Confirmation						◆		
70111	Gabapentin, Urine (Confirmation)							◆	
74165	Gaba, Urine Confirmation						◆		
70153	Heroin, Urine (Confirmation)							◆	
74028	Heroin, Urine Confirmation						◆		
74170	Hydrocodone, Urine Confirmation						◆		
74180	Imipramine, Urine Confirmation						◆		
74190	Ketamine, Urine Confirmation						◆		
2570	Lipoprotein (a)							◆	
92203	Lipoprotein (a)						◆		
70113	Meperidine, Urine (Confirmation)							◆	
74119	Meperidine, Urine Confirmation						◆		
70115	Methadone, Urine (Confirmation)							◆	
74200	Methadone, Urine Confirmation						◆		
70154	Methylphenidate, Urine (Confirmation)							◆	
74210	Methylphenidate, Urine Confirmation						◆		
70125	Naloxone, Urine (Confirmation)							◆	
74044	Naloxone, Urine Confirmation						◆		
74220	Naltrexone, Urine Confirmation						◆		
91352	Nicotine, Urine							◆	
74230	Nicotine, Urine Confirmation						◆		
91610	NMR Lipoprofile							◆	
92186	NMR Lipoprofile						◆		
70118	Opiates, Urine (Confirmation)							◆	
74240	Opiates, Urine Confirmation						◆		
70420	Oxycodone, Urine (Confirmation)							◆	
74250	Oxycodone, Urine Confirmation						◆		
70126	Phencyclidine, Urine (Confirmation)							◆	
74027	Phencyclidine, Urine Confirmation						◆		
74054	Propoxyphene, Urine Confirmation						◆		
74260	Suboxone, Urine Confirmation						◆		
70157	Tapentadol, Urine (Confirmation)							◆	
74290	Tapentadol, Urine Confirmation						◆		
70135	Tramadol, Urine (Confirmation)							◆	
74270	Tramadol, Urine Confirmation						◆		

70077 Amphetamine, Urine (Confirmation)

DT

74020 Amphetamine, Urine Confirmation

NT

Specimen:	
Collect:	Random Urine in Sterile Specimen Container
Submit:	10 mL (Min:3 mL) Random Urine. Submit Refrigerated.
Special Handling:	Avoid Repeated Freeze/Thaw Cycles
Rejection Criteria:	Urine stabilized with additives or preservatives
Stability:	Ambient: 1 Week(s); Refrigerated: 11 Day(s); Frozen: 2 Week(s); Incubated: Unacceptable
Methodology:	Liquid Chromatography Mass Spectrometry
Performed:	Mon-Fri
Reported:	1-6 Day(s)
CPT Codes:	80324
Interpretive Data:	Please see report for interpretive data.
Components:	74059 - AMPHETAMINE 74113 - MDEA 74064 - METHAMPHET 74114 - PHENTERMINE

70128 Antidepressant (SSRI), Urine (Confirmation)

DT

74030 Antidepressant (SSRI), Urine Confirmation

NT

Specimen:	
Collect:	Random Urine in Sterile Specimen Container
Submit:	10 mL (Min:3 mL) Random Urine. Submit Refrigerated.
Special Handling:	Avoid Repeated Freeze/Thaw Cycles
Rejection Criteria:	Urine stabilized with additives or preservatives
Stability:	Ambient: 1 Week(s); Refrigerated: 2 Week(s); Frozen: 2 Week(s); Incubated: Unacceptable
Methodology:	Liquid Chromatography Mass Spectrometry
Performed:	Mon-Fri
Reported:	1-6 Day(s)
CPT Codes:	80336
Interpretive Data:	Please see report for interpretive data.
Components:	74083 - FLUOXETINE 74093 - SERTRALINE 74076 - CITALOPRAM 74092 - PAROXETINE 74088 - NORFLUOXETINE 74082 - DULOXETINE 74087 - NDESMETHYLCIT

74040 Antidepressant (SSRI 2), Urine Confirmation
NT

Specimen:					
Collect:	Random Urine in Sterile Specimen Container				
Submit:	10 mL (Min:3 mL) Random Urine. Submit Refrigerated.				
Special Handling:	Avoid Repeated Freeze/Thaw Cycles				
Rejection Criteria:	Urine stabilized with additives or preservatives				
Stability:	Ambient: 1 Week(s); Refrigerated: 2 Week(s); Frozen: 2 Week(s); Incubated: Unacceptable				
Methodology:	Liquid Chromatography Mass Spectrometry				
Performed:	Mon-Fri				
Reported:	1-6 Day(s)				
CPT Codes:	80336				
Interpretive Data:	Please see report for interpretive data.				
Components:	<table style="width: 100%; border: none;"> <tr> <td style="width: 50%;">74074 - BUPROPION</td> <td style="width: 50%;">74084 - OHBUPROPION</td> </tr> <tr> <td>74091 - ODESMETVENLAF</td> <td>74096 - VENLAFAXINE</td> </tr> </table>	74074 - BUPROPION	74084 - OHBUPROPION	74091 - ODESMETVENLAF	74096 - VENLAFAXINE
74074 - BUPROPION	74084 - OHBUPROPION				
74091 - ODESMETVENLAF	74096 - VENLAFAXINE				

74050 Antidepressant (TCA), Urine Confirmation
NT

Specimen:									
Collect:	Random Urine in Sterile Specimen Container								
Submit:	10 mL (Min:3 mL) Random Urine. Submit Refrigerated.								
Special Handling:	Avoid Repeated Freeze/Thaw Cycles								
Rejection Criteria:	Urine stabilized with additives or preservatives								
Stability:	Ambient: 1 Week(s); Refrigerated: 2 Week(s); Frozen: 2 Week(s); Incubated: Unacceptable								
Methodology:	Liquid Chromatography Mass Spectrometry								
Performed:	Mon-Fri								
Reported:	1-6 Day(s)								
CPT Codes:	80337								
Interpretive Data:	Please see report for interpretive data.								
Components:	<table style="width: 100%; border: none;"> <tr> <td style="width: 50%;">74073 - AMITRIPTYLINE</td> <td style="width: 50%;">74077 - CLOMIPRAMINE</td> </tr> <tr> <td>74078 - DESIPRAMINE</td> <td>74079 - DESMETDOXEPIN</td> </tr> <tr> <td>74081 - DOXEPIN</td> <td>74086 - IMIPRAMINE</td> </tr> <tr> <td>74089 - NORTRIPTYLINE</td> <td>74094 - TRIMIPRAMINE</td> </tr> </table>	74073 - AMITRIPTYLINE	74077 - CLOMIPRAMINE	74078 - DESIPRAMINE	74079 - DESMETDOXEPIN	74081 - DOXEPIN	74086 - IMIPRAMINE	74089 - NORTRIPTYLINE	74094 - TRIMIPRAMINE
74073 - AMITRIPTYLINE	74077 - CLOMIPRAMINE								
74078 - DESIPRAMINE	74079 - DESMETDOXEPIN								
74081 - DOXEPIN	74086 - IMIPRAMINE								
74089 - NORTRIPTYLINE	74094 - TRIMIPRAMINE								

70080 Barbiturates, Urine (Confirmation)

DT

74060 Barbiturates, Urine Confirmation

NT

Specimen:							
Collect:	Random Urine in Sterile Specimen Container						
Submit:	10 mL (Min:3 mL) Random Urine. Submit Refrigerated.						
Special Handling:	Avoid Repeated Freeze/Thaw Cycles						
Rejection Criteria:	Urine stabilized with additives or preservatives						
Stability:	Ambient: 1 Week(s); Refrigerated: 1 Month(s); Frozen: 24 Month(s); Incubated: Unacceptable						
Methodology:	Liquid Chromatography Mass Spectrometry						
Performed:	Mon-Fri						
Reported:	1-6 Day(s)						
CPT Codes:	80345						
Interpretive Data:	Please see report for interpretive data.						
Components:	<table border="0"> <tr> <td>74004 - BUTALBITAL</td> <td>74007 - SECOBARBITAL</td> </tr> <tr> <td>74006 - PHENOBARBITAL</td> <td>74132 - PENTOBARBITAL</td> </tr> <tr> <td>74133 - AMOBARBITAL</td> <td></td> </tr> </table>	74004 - BUTALBITAL	74007 - SECOBARBITAL	74006 - PHENOBARBITAL	74132 - PENTOBARBITAL	74133 - AMOBARBITAL	
74004 - BUTALBITAL	74007 - SECOBARBITAL						
74006 - PHENOBARBITAL	74132 - PENTOBARBITAL						
74133 - AMOBARBITAL							

70089 Benzodiazepines, Urine (Confirmation)

DT

74070 Benzodiazepine, Urine Confirmation

NT

Specimen:											
Collect:	Random Urine in Sterile Specimen Container										
Submit:	10 mL (Min:3 mL) Random Urine. Submit Refrigerated.										
Special Handling:	Avoid Repeated Freeze/Thaw Cycles										
Rejection Criteria:	Urine stabilized with additives or preservatives										
Stability:	Ambient: 1 Week(s); Refrigerated: 1 Week(s); Frozen: 1 Week(s); Incubated: Unacceptable										
Methodology:	Liquid Chromatography Mass Spectrometry										
Performed:	Mon-Fri										
Reported:	1-6 Day(s)										
CPT Codes:	80346										
Interpretive Data:	Please see report for interpretive data.										
Components:	<table border="0"> <tr> <td>74008 - 7AMCLONAZEPAM</td> <td>74009 - AOHALPRAZOLAM</td> </tr> <tr> <td>74012 - LORAZEPAM</td> <td>74013 - NORDIAZEPAM</td> </tr> <tr> <td>74014 - OXAZEPAM</td> <td>74016 - TEMAZEPAM</td> </tr> <tr> <td>74116 - ALPRAZOLAM</td> <td>74117 - DIAZEPAM</td> </tr> <tr> <td>74118 - MIDAZOLAM</td> <td>74011 - CLONAZEPAM</td> </tr> </table>	74008 - 7AMCLONAZEPAM	74009 - AOHALPRAZOLAM	74012 - LORAZEPAM	74013 - NORDIAZEPAM	74014 - OXAZEPAM	74016 - TEMAZEPAM	74116 - ALPRAZOLAM	74117 - DIAZEPAM	74118 - MIDAZOLAM	74011 - CLONAZEPAM
74008 - 7AMCLONAZEPAM	74009 - AOHALPRAZOLAM										
74012 - LORAZEPAM	74013 - NORDIAZEPAM										
74014 - OXAZEPAM	74016 - TEMAZEPAM										
74116 - ALPRAZOLAM	74117 - DIAZEPAM										
74118 - MIDAZOLAM	74011 - CLONAZEPAM										

70086 Buprenorphine, Urine (Confirmation)
DT
74080 Buprenorphine, Urine Confirmation
NT

Specimen:	
Collect:	Random Urine in Sterile Specimen Container
Submit:	10 mL (Min:3 mL) Random Urine. Submit Refrigerated.
Special Handling:	Avoid Repeated Freeze/Thaw Cycles
Rejection Criteria:	Urine stabilized with additives or preservatives
Stability:	Ambient: 1 Week(s); Refrigerated: 1 Month(s); Frozen: 24 Month(s); Incubated: Unacceptable
Methodology:	Liquid Chromatography Mass Spectrometry
Performed:	Mon-Fri
Reported:	1-6 Day(s)
CPT Codes:	80348
Interpretive Data:	Please see report for interpretive data.
Components:	74031 - BUPRENORPHINE 74047 - NORBUPRENORPH

70133 Cannabinoids, Urine (Confirmation)
DT
74017 Cannabinoids, Urine Confirmation
NT

Specimen:	
Collect:	Random Urine in Sterile Specimen Container
Submit:	10 mL (Min:3 mL) Random Urine. Submit Refrigerated.
Special Handling:	Avoid Repeated Freeze/Thaw Cycles
Rejection Criteria:	Urine stabilized with additives or preservatives
Stability:	Ambient: 1 Week(s); Refrigerated: 1 Month(s); Frozen: 24 Month(s); Incubated: Unacceptable
Methodology:	Liquid Chromatography Mass Spectrometry
Performed:	Mon-Fri
Reported:	1-6 Day(s)
CPT Codes:	80349
Interpretive Data:	Please see report for interpretive data.

2759 Cardiovascular Disease Risk Panel + NMR, Lipoprofile CC

Specimen:		
Collect:	One LipoTube One SST	
Submit:	2 mL (Min:2 mL) Plasma in LipoTube. Submit Refrigerated. 1 mL (Min:0.5 mL) Serum. Submit Refrigerated. Submit in a Standard Transport Tube.	
Special Handling:	See Components	
Rejection Criteria:	See Components	
Stability:	Ambient: 6 Hour(s); Refrigerated: 6 Day(s); Frozen: Unacceptable; Incubated: Unacceptable	
Methodology:	See Individual Components	
Performed:	Mon-Fri	
Reported:	4-6 Day(s)	
CPT Codes:	80061 83704	83695 83721
Components:	2032 - LDL DIRECT 91610 - NMR LipoProfile 92186 - NMR LipoProfile	2570 - LIOPROTEIN (a) 92203 - LIOPROTEIN (a)

Please take note of changes to components.

Component Changes:

Remove: 2570 - LIOPROTEIN (a); 91610 - NMR LipoProfile

Add: 92203 - LIOPROTEIN (a); 92186 - NMR LipoProfile

70101 Carisoprodol, Urine (Confirmation) DT
74090 Carisoprodol, Urine Confirmation NT

Specimen:		
Collect:	Random Urine in Sterile Specimen Container	
Submit:	10 mL (Min:3 mL) Random Urine. Submit Refrigerated.	
Special Handling:	Avoid Repeated Freeze/Thaw Cycles	
Rejection Criteria:	Urine stabilized with additives or preservatives	
Stability:	Ambient: 3 Week(s); Refrigerated: 3 Week(s); Frozen: 3 Week(s); Incubated: Unacceptable	
Methodology:	Liquid Chromatography Mass Spectrometry	
Performed:	Mon-Fri	
Reported:	1-6 Day(s)	
CPT Codes:	80369	
Interpretive Data:	Please see report for interpretive data.	
Components:	74021 - CARISOPRODOL	74023 - MEPROBAMATE

70100 Cocaine, Urine (Confirmation)
DT
74062 Cocaine, Urine Confirmation
NT

Specimen:	
Collect:	Random Urine in Sterile Specimen Container
Submit:	10 mL (Min:3 mL) Random Urine. Submit Refrigerated.
Special Handling:	Avoid Repeated Freeze/Thaw Cycles
Rejection Criteria:	Urine stabilized with additives or preservatives
Stability:	Ambient: 1 Week(s); Refrigerated: 1 Month(s); Frozen: 24 Month(s); Incubated: Unacceptable
Methodology:	Liquid Chromatography Mass Spectrometry
Performed:	Mon-Fri
Reported:	1-6 Day(s)
CPT Codes:	80353
Interpretive Data:	Please see report for interpretive data.

74100 Codeine + Morphine, Urine Confirmation
NT

Specimen:	
Collect:	Random Urine in Sterile Specimen Container
Submit:	10 mL (Min:3 mL) Random Urine. Submit Refrigerated.
Special Handling:	Avoid Repeated Freeze/Thaw Cycles
Rejection Criteria:	Urine stabilized with additives or preservatives
Stability:	Ambient: 1 Week(s); Refrigerated: 1 Month(s); Frozen: 24 Month(s); Incubated: Unacceptable
Methodology:	Liquid Chromatography Mass Spectrometry
Performed:	Mon-Fri
Reported:	1-6 Day(s)
CPT Codes:	80361
Interpretive Data:	Please see report for interpretive data.
Components:	74032 - CODEINE 74043 - MORPHINE

70104 Cyclobenzaprine, Urine (Confirmation)
DT
74022 Cyclobenzaprine, Urine Confirmation
NT

Specimen:	
Collect:	Random Urine in Sterile Specimen Container
Submit:	10 mL (Min:3 mL) Random Urine. Submit Refrigerated.
Special Handling:	Avoid Repeated Freeze/Thaw Cycles
Rejection Criteria:	Urine stabilized with additives or preservatives
Stability:	Ambient: 1 Month(s); Refrigerated: 1 Month(s); Frozen: 1 Month(s); Incubated: Unacceptable
Methodology:	Liquid Chromatography Mass Spectrometry
Performed:	Mon-Fri
Reported:	1-6 Day(s)
CPT Codes:	80335
Interpretive Data:	Please see report for interpretive data.

74110 Dextromethorphan, Urine Confirmation
NT

Specimen:	
Collect:	Random Urine in Sterile Specimen Container
Submit:	10 mL (Min:3 mL) Random Urine. Submit Refrigerated.
Special Handling:	Avoid Repeated Freeze/Thaw Cycles
Rejection Criteria:	Urine stabilized with additives or preservatives
Stability:	Ambient: 2 Hour(s); Refrigerated: 3 Day(s); Frozen: 12 Month(s); Incubated: Unacceptable
Methodology:	Liquid Chromatography Mass Spectrometry
Performed:	Mon-Fri
Reported:	1-6 Day(s)
CPT Codes:	80362
Interpretive Data:	Please see report for interpretive data.
Components:	74003 - DEXMETHORPHAN 74038 - LEVORPHANOL

74120 Doxepin, Urine Confirmation
NT

Specimen:	
Collect:	Random Urine in Sterile Specimen Container
Submit:	10 mL (Min:3 mL) Random Urine. Submit Refrigerated.
Special Handling:	Avoid Repeated Freeze/Thaw Cycles
Rejection Criteria:	Urine stabilized with additives or preservatives
Stability:	Ambient: 1 Week(s); Refrigerated: 2 Week(s); Frozen: 2 Week(s); Incubated: Unacceptable
Methodology:	Liquid Chromatography Mass Spectrometry
Performed:	Mon-Fri
Reported:	1-6 Day(s)
CPT Codes:	80335
Interpretive Data:	Please see report for interpretive data.
Components:	74081 - DOXEPIN 74079 - DESMETDOXEPIN

6955 Drug Screen, Urine Medical
DT

73010 Drug Screen, Urine Medication Management
NT/AOE

Specimen:																	
Collect:	Random Urine in Sterile Specimen Container																
Submit:	10 mL (Min:3 mL) Random Urine. Submit Refrigerated.																
Special Handling:	<p>Avoid Repeated Freeze/Thaw Cycles Recommended pH is 3.0-11.0</p> <p>IF POSITIVE CONFIRMATION WILL BE PERFORMED. Additional charges will apply. ALCOHOL: Ethyl Glucuronide, Ethyl Sulfate AMPHETAMINES: Amphetamine, Methamphetamine, MDEA-Eve, Phentermine BARBITURATES: Amobarbital, Butalbital, Pentobarbital, Phenobarbital, Secobarbital BENZODIAZEPINES: 7-Aminoclonazepam, Alpha-hydroxyalprazolam, Alprazolam, Clonazepam, Lorazepam, Midazolam, Nordiazepam, Oxazepam BUPRENORPHINE: Buprenorphine, Norbuprenorphine CANNABINOIDS: 11-Nor-9-carboxy-THC COCAINE: Benzoylecgonine ECSTASY: Methylenedioxyamphetamine (MDA), Methylenedioxymethamphetamine (MDMA) FENTANYL: Acetylfentanyl, Acetylnorfentanyl, Fentanyl, Norfentanyl METHADONE: EDDP, Methadone OPIATES: 6-Acetylmorphine, Codeine, Desmethyltramadol, Hydrocodone, Hydromorphone, Norhydrocodone, Noroxycodone, Morphine, Tapentadol-o-Sulf, Tramadol, Zolpidem OXYCODONE: Noroxycodone, Oxycodone, Oxycodone, Oxycodone PHENCYCLIDINE: Phencyclidine PROPOXYPHENE: Propoxyphene</p>																
Rejection Criteria:	Urine stabilized with additives or preservatives																
Stability:	Ambient: Unacceptable; Refrigerated: 3 Day(s); Frozen: 1 Week(s); Incubated: Unacceptable																
Methodology:	Homogeneous Enzyme Immunoassay; Liquid Chromatography Mass Spectrometry																
Performed:	Mon-Fri																
Reported:	1-6 Day(s)																
CPT Codes:	80307 82570																
Interpretive Data:	Please see report for interpretive data.																
Components:	<table border="0"> <tr> <td>7001 - AMPHETAMINES</td> <td>7002 - BARBITURATES</td> </tr> <tr> <td>7003 - BENZODIAZEPINES</td> <td>7004 - CANNABINOIDS</td> </tr> <tr> <td>7005 - COCAINE</td> <td>7006 - ECSTASY</td> </tr> <tr> <td>7007 - METHADONE</td> <td>7008 - OPIATES</td> </tr> <tr> <td>7009 - PHENCYCLIDINE</td> <td>7011 - OXYCODONE</td> </tr> <tr> <td>7012 - CREATININE, URINE</td> <td>7013 - ALCOHOL</td> </tr> <tr> <td>73031 - BUPRENORPHINE</td> <td>73034 - FENTANYL</td> </tr> <tr> <td>73054 - PROPOXYPHENE</td> <td></td> </tr> </table>	7001 - AMPHETAMINES	7002 - BARBITURATES	7003 - BENZODIAZEPINES	7004 - CANNABINOIDS	7005 - COCAINE	7006 - ECSTASY	7007 - METHADONE	7008 - OPIATES	7009 - PHENCYCLIDINE	7011 - OXYCODONE	7012 - CREATININE, URINE	7013 - ALCOHOL	73031 - BUPRENORPHINE	73034 - FENTANYL	73054 - PROPOXYPHENE	
7001 - AMPHETAMINES	7002 - BARBITURATES																
7003 - BENZODIAZEPINES	7004 - CANNABINOIDS																
7005 - COCAINE	7006 - ECSTASY																
7007 - METHADONE	7008 - OPIATES																
7009 - PHENCYCLIDINE	7011 - OXYCODONE																
7012 - CREATININE, URINE	7013 - ALCOHOL																
73031 - BUPRENORPHINE	73034 - FENTANYL																
73054 - PROPOXYPHENE																	

Please take note of ask on order entry questions

1136 List prescribed medications

1164 Do you want all prescribed medications confirmed?

73020 Drug Screen, Urine Medication Management w/o THC and Alcohol NT/AOE

Specimen:															
Collect:	Random Urine in Sterile Specimen Container														
Submit:	10 mL (Min:3 mL) Random Urine. Submit Refrigerated.														
Special Handling:	Avoid Repeated Freeze/Thaw Cycles Recommended pH is 3.0-11.0 IF POSITIVE THEN CONFIRMATION WILL BE PERFORMED. Additional charges will apply. AMPHETAMINES: Amphetamine, Methamphetamine, MDEA-Eve, Phentermine BARBITURATES: Amobarbital, Butalbital, Pentobarbital, Phenobarbital, Secobarbital BENZODIAZEPINES: 7-Aminoclonazepam, Alpha-hydroxyalprazolam, Alprazolam, Clonazepam, Lorazepam, Midazolam, Nordiazepam, Oxazepam BUPRENORPHINE: Buprenorphine, Norbuprenorphine COCAINE: Benzoylcegonine ECSTASY: Methylenedioxyamphetamine (MDA), Methylenedioxymethamphetamine (MDMA) FENTANYL: Acetylfentanyl, Acetylnorfentanyl, Fentanyl, Norfentanyl METHADONE: EDDP, Methadone OPIATES: 6-Acetylmorphine, Codeine, Hydrocodone, Hydromorphone, Norhydrocodone, Noroxycodone, Morphine, Tapentadol-o-Sulf, Tramadol, Zolpidem OXYCODONE: Noroxycodone, Oxycodone, Oxycodone, Oxycodone PHENCYCLIDINE: Phencyclidine PROPOXYPHENE: Propoxyphene														
Rejection Criteria:	Urine stabilized with additives or preservatives														
Stability:	Ambient: Unacceptable; Refrigerated: 3 Day(s); Frozen: 1 Week(s); Incubated: Unacceptable														
Methodology:	Homogeneous Enzyme Immunoassay; Liquid Chromatography Mass Spectrometry														
Performed:	Mon-Fri														
Reported:	1-6 Day(s)														
CPT Codes:	80307 82570														
Interpretive Data:	Please see report for interpretive data.														
Components:	<table border="0"> <tr> <td>7001 - AMPHETAMINES</td> <td>7002 - BARBITURATES</td> </tr> <tr> <td>7003 - BENZODIAZEPINES</td> <td>7005 - COCAINE</td> </tr> <tr> <td>7006 - ECSTASY</td> <td>7007 - METHADONE</td> </tr> <tr> <td>7008 - OPIATES</td> <td>7009 - PHENCYCLIDINE</td> </tr> <tr> <td>7011 - OXYCODONE</td> <td>7012 - CREATININE, URINE</td> </tr> <tr> <td>73031 - BUPRENORPHINE</td> <td>73034 - FENTANYL</td> </tr> <tr> <td>73054 - PROPOXYPHENE</td> <td></td> </tr> </table>	7001 - AMPHETAMINES	7002 - BARBITURATES	7003 - BENZODIAZEPINES	7005 - COCAINE	7006 - ECSTASY	7007 - METHADONE	7008 - OPIATES	7009 - PHENCYCLIDINE	7011 - OXYCODONE	7012 - CREATININE, URINE	73031 - BUPRENORPHINE	73034 - FENTANYL	73054 - PROPOXYPHENE	
7001 - AMPHETAMINES	7002 - BARBITURATES														
7003 - BENZODIAZEPINES	7005 - COCAINE														
7006 - ECSTASY	7007 - METHADONE														
7008 - OPIATES	7009 - PHENCYCLIDINE														
7011 - OXYCODONE	7012 - CREATININE, URINE														
73031 - BUPRENORPHINE	73034 - FENTANYL														
73054 - PROPOXYPHENE															

Please take note of ask on order entry questions.

1136 List prescribed medications

1164 Do you want all prescribed medications confirmed?

7000 Drug Screen, Urine Pain Management
DT
73000 Drug Screen, Urine Pain Management
NT

Specimen:	
Collect:	Random Urine in Sterile Specimen Container
Submit:	10 mL (Min:3 mL) Random Urine. Submit Refrigerated.
Special Handling:	Avoid Repeated Freeze/Thaw Cycles Recommended pH is 3.0-11.0 IF POSITIVE THEN CONFIRMATION WILL BE PERFORMED. Additional charges will apply. ALCOHOL: Ethyl Glucuronide, Ethyl Sulfate AMPHETAMINES: Amphetamine, Methamphetamine, MDEA-Eve, Phentermine BARBITURATES: Amobarbital, Butalbital, Pentobarbital, Phenobarbital, Secobarbital BENZODIAZEPINES: 7-Aminoclonazepam, Alpha-hydroxyalprazolam, Alprazolam, Clonazepam, Lorazepam, Midazolam, Nordiazepam, Oxazepam CANNABINOIDS: 11-Nor-9-carboxy-THC COCAINE: Benzoylecgonine ECSTASY: Methylenedioxyamphetamine (MDA), Methylenedioxymethamphetamine (MDMA) METHADONE: EDDP, Methadone OPIATES: 6-Acetylmorphine, Codeine, Desmethyltramadol, Hydrocodone, Hydromorphone, Norhydrocodone, Noroxymorphone, Morphine, Tapentadol-o-Sulf, Tramadol, Zolpidem OXYCODONE: Noroxycodone, Oxycodone, Oxymorphone PHENCYCLIDINE: Phencyclidine
Rejection Criteria:	Urine stabilized with additives or preservatives
Stability:	Ambient: 2 Hour(s); Refrigerated: 3 Day(s); Frozen: 12 Month(s); Incubated: Unacceptable
Methodology:	Homogeneous Enzyme Immunoassay
Performed:	Sunday, Mon-Fri
Reported:	1-2 Day(s)
CPT Codes:	80307 82570
Interpretive Data:	Please see report for interpretive data.
Components:	7001 - AMPHETAMINES 7003 - BENZODIAZEPINES 7005 - COCAINE 7007 - METHADONE 7009 - PHENCYCLIDINE 7012 - CREATININE, URINE 7002 - BARBITURATES 7004 - CANNABINOIDS 7006 - ECSTASY 7008 - OPIATES 7011 - OXYCODONE 7013 - ALCOHOL

70137 Ecstasy, Urine (Confirmation)
DT
74130 Ecstasy, Urine Confirmation
NT

Specimen:	
Collect:	Random Urine in Sterile Specimen Container
Submit:	10 mL (Min:3 mL) Random Urine. Submit Refrigerated.
Special Handling:	Avoid Repeated Freeze/Thaw Cycles
Rejection Criteria:	Urine stabilized with additives or preservatives
Stability:	Ambient: 1 Week(s); Refrigerated: 1 Month(s); Frozen: 24 Month(s); Incubated: Unacceptable
Methodology:	Liquid Chromatography Mass Spectrometry
Performed:	Mon-Fri
Reported:	1-6 Day(s)
CPT Codes:	80359
Interpretive Data:	Please see report for interpretive data.
Components:	74066 - MDAMPHETAMINE 74067 - MDMETHAMPHET

74140 Embeda, Urine Confirmation
NT

Specimen:	
Collect:	Random Urine in Sterile Specimen Container
Submit:	10 mL (Min:3 mL) Random Urine. Submit Refrigerated.
Special Handling:	Avoid Repeated Freeze/Thaw Cycles
Rejection Criteria:	Urine stabilized with additives or preservatives
Stability:	Ambient: 5 Day(s); Refrigerated: 1 Week(s); Frozen: 1 Month(s); Incubated: Unacceptable
Methodology:	Liquid Chromatography Mass Spectrometry
Performed:	Mon-Fri
Reported:	1-6 Day(s)
CPT Codes:	80361
Interpretive Data:	Please see report for interpretive data.
Components:	74043 - MORPHINE 74029 - 6BNALTREXOL 74046 - NALTREXONE

70111 Gabapentin, Urine (Confirmation)
DT
74165 Gaba, Urine Confirmation
NT

Specimen:	
Collect:	Random Urine in Sterile Specimen Container
Submit:	10 mL (Min:3 mL) Random Urine. Submit Refrigerated.
Special Handling:	Avoid Repeated Freeze/Thaw Cycles
Rejection Criteria:	Urine stabilized with additives or preservatives
Stability:	Ambient: 1 Month(s); Refrigerated: 1 Month(s); Frozen: 1 Month(s); Incubated: Unacceptable
Methodology:	Liquid Chromatography Mass Spectrometry
Performed:	Mon-Fri
Reported:	1-6 Day(s)
Interpretive Data:	Please see report for interpretive data.
Components:	74001 - GABAPENTIN 74002 - PREGABALIN

70153 Heroin, Urine (Confirmation)
DT
74028 Heroin, Urine Confirmation
NT

Specimen:	
Collect:	Random Urine in Sterile Specimen Container
Submit:	10 mL (Min:3 mL) Random Urine. Submit Refrigerated.
Special Handling:	Avoid Repeated Freeze/Thaw Cycles
Rejection Criteria:	Urine stabilized with additives or preservatives
Stability:	Ambient: 1 Week(s); Refrigerated: 1 Month(s); Frozen: 24 Month(s); Incubated: Unacceptable
Methodology:	Liquid Chromatography Mass Spectrometry
Performed:	Mon-Fri
Reported:	1-6 Day(s)
CPT Codes:	80356
Interpretive Data:	<p>General Reference Range : <10 ng/mL Values listed in the reference/normal range indicate a laboratory reporting cutoff.</p> <p>Result below the cutoff value (negative) may indicate non-use, poor timing of specimen collection, poor absorption, specimen dilution, specimen adulteration, or testing limitation. Quantitative results are reported only when equal to or above cutoff value.</p> <p>For medical purpose only. Methodology: Quantitative Liquid Chromatography Tandem Mass Spectrometry. Test developed and validated by Interpath Laboratory. Questions can be directed to laboratory.</p> <p>Laboratory Director: Jane J Yin MD, PhD</p>

74170 Hydrocodone, Urine Confirmation
NT

Specimen:	
Collect:	Random Urine in Sterile Specimen Container
Submit:	10 mL (Min:3 mL) Random Urine. Submit Refrigerated.
Special Handling:	Avoid Repeated Freeze/Thaw Cycles
Rejection Criteria:	Urine stabilized with additives or preservatives
Stability:	Ambient: 1 Week(s); Refrigerated: 1 Month(s); Frozen: 24 Month(s); Incubated: Unacceptable
Methodology:	Liquid Chromatography Mass Spectrometry
Performed:	Mon-Fri
Reported:	1-6 Day(s)
CPT Codes:	80361
Interpretive Data:	Please see report for interpretive data.
Components:	74036 - HYDROCODONE 74049 - NORHYDROCODON

74180 Imipramine, Urine Confirmation
NT

Specimen:	
Collect:	Random Urine in Sterile Specimen Container
Submit:	10 mL (Min:3 mL) Random Urine. Submit Refrigerated.
Special Handling:	Avoid Repeated Freeze/Thaw Cycles
Rejection Criteria:	Urine stabilized with additives or preservatives
Stability:	Ambient: 1 Week(s); Refrigerated: 2 Week(s); Frozen: 2 Week(s); Incubated: Unacceptable
Methodology:	Liquid Chromatography Mass Spectrometry
Performed:	Mon-Fri
Reported:	1-6 Day(s)
CPT Codes:	80335
Interpretive Data:	Please see report for interpretive data.
Components:	74078 - DESIPRAMINE 74086 - IMIPRAMINE 74107 - IMI+DESIPRAM

74190 Ketamine, Urine Confirmation
NT

Specimen:	
Collect:	Random Urine in Sterile Specimen Container
Submit:	10 mL (Min:3 mL) Random Urine. Submit Refrigerated.
Special Handling:	Avoid Repeated Freeze/Thaw Cycles
Rejection Criteria:	Urine stabilized with additives or preservatives
Stability:	Ambient: 2 Day(s); Refrigerated: 2 Week(s); Frozen: 2 Week(s); Incubated: Unacceptable
Methodology:	Liquid Chromatography Mass Spectrometry
Performed:	Mon-Fri
Reported:	1-6 Day(s)
CPT Codes:	80357
Interpretive Data:	Please see report for interpretive data.
Components:	74024 - KETAMINE 74026 - NORKETAMINE

2570 Lipoprotein (a)
DT
92203 Lipoprotein (a)
NT

Specimen:	
Collect:	One SST Also Acceptable One Green Top (Li Heparin) One Lavender (EDTA) One Pink Top (EDTA)
Submit:	1 mL (Min:0.5 mL) Serum. Submit Frozen. Submit in a Standard Transport Tube. Also Acceptable 1 mL (Min:0.5 mL) Plasma. Submit Refrigerated. Submit in a Standard Transport Tube.
Stability:	Ambient: 8 Hour(s); Refrigerated: 2 Day(s); Frozen: 3 Month(s); Incubated: Unacceptable
Methodology:	Quantitative Immunoturbidimetry
Performed:	Sun-Sat
Reported:	2-4 Day(s)
CPT Codes:	83695

70113 Meperidine, Urine (Confirmation)
DT
74119 Meperidine, Urine Confirmation
NT

Specimen:	
Collect:	Random Urine in Sterile Specimen Container
Submit:	10 mL (Min:3 mL) Random Urine. Submit Refrigerated.
Special Handling:	Avoid Repeated Freeze/Thaw Cycles
Rejection Criteria:	Urine stabilized with additives or preservatives
Stability:	Ambient: 1 Week(s); Refrigerated: 1 Month(s); Frozen: 24 Month(s); Incubated: Unacceptable
Methodology:	Liquid Chromatography Mass Spectrometry
Performed:	Mon-Fri
Reported:	1-6 Day(s)
CPT Codes:	80361
Interpretive Data:	Please see report for interpretive data.

70115 Methadone, Urine (Confirmation)
DT
74200 Methadone, Urine Confirmation
NT

Specimen:	
Collect:	Random Urine in Sterile Specimen Container
Submit:	10 mL (Min:3 mL) Random Urine. Submit Refrigerated.
Special Handling:	Avoid Repeated Freeze/Thaw Cycles
Rejection Criteria:	Urine stabilized with additives or preservatives
Stability:	Ambient: 1 Week(s); Refrigerated: 1 Month(s); Frozen: 24 Month(s); Incubated: Unacceptable
Methodology:	Liquid Chromatography Mass Spectrometry
Performed:	Mon-Fri
Reported:	1-6 Day(s)
CPT Codes:	80358
Interpretive Data:	Please see report for interpretive data.
Components:	74033 - METHADONEMET 74041 - METHADONE

70154 Methylphenidate, Urine (Confirmation)
DT
74210 Methylphenidate, Urine Confirmation
NT

Specimen:	
Collect:	Random Urine in Sterile Specimen Container
Submit:	10 mL (Min:3 mL) Random Urine. Submit Refrigerated.
Special Handling:	Avoid Repeated Freeze/Thaw Cycles
Rejection Criteria:	Urine stabilized with additives or preservatives
Stability:	Ambient: Unacceptable; Refrigerated: 3 Week(s); Frozen: 3 Month(s); Incubated: Unacceptable
Methodology:	Liquid Chromatography Mass Spectrometry
Performed:	Mon-Fri
Reported:	1-6 Day(s)
CPT Codes:	80360
Interpretive Data:	Please see report for interpretive data.
Components:	74068 - MPHENIDATE 74072 - RITALINICACID

70125 Naloxone, Urine (Confirmation)
DT
74044 Naloxone, Urine Confirmation
NT

Specimen:	
Collect:	Random Urine in Sterile Specimen Container
Submit:	10 mL (Min:3 mL) Random Urine. Submit Refrigerated.
Special Handling:	Avoid Repeated Freeze/Thaw Cycles
Rejection Criteria:	Urine stabilized with additives or preservatives
Stability:	Ambient: 1 Week(s); Refrigerated: 1 Month(s); Frozen: 24 Month(s); Incubated: Unacceptable
Methodology:	Liquid Chromatography Mass Spectrometry
Performed:	Mon-Fri
Reported:	1-6 Day(s)
CPT Codes:	80362
Interpretive Data:	Please see report for interpretive data.

74220 Naltrexone, Urine Confirmation
NT

Specimen:	
Collect:	Random Urine in Sterile Specimen Container
Submit:	10 mL (Min:3 mL) Random Urine. Submit Refrigerated.
Special Handling:	Avoid Repeated Freeze/Thaw Cycles
Rejection Criteria:	Urine stabilized with additives or preservatives
Stability:	Ambient: 5 Day(s); Refrigerated: 1 Week(s); Frozen: 1 Month(s); Incubated: Unacceptable
Methodology:	Liquid Chromatography Mass Spectrometry
Performed:	Mon-Fri
Reported:	1-6 Day(s)
CPT Codes:	80362
Interpretive Data:	Please see report for interpretive data.
Components:	74029 - 6BNALTREXOL 74046 - NALTREXONE

91352 Nicotine, Urine
DT
74230 Nicotine, Urine Confirmation
NT

Specimen:							
Collect:	Random Urine in Sterile Specimen Container						
Submit:	10 mL (Min:3 mL) Random Urine. Submit Refrigerated.						
Special Handling:	Avoid Repeated Freeze/Thaw Cycles						
Rejection Criteria:	Urine stabilized with additives or preservatives						
Stability:	Ambient: 10 Day(s); Refrigerated: 10 Day(s); Frozen: 8 Month(s); Incubated: Unacceptable						
Methodology:	Liquid Chromatography Mass Spectrometry						
Performed:	Mon-Fri						
Reported:	1-6 Day(s)						
CPT Codes:	80323						
Interpretive Data:	Please see report for interpretive data.						
Components:	<table style="width: 100%; border: none;"> <tr> <td style="width: 50%;">74069 - NICOTINE</td> <td style="width: 50%;">74063 - COTININE</td> </tr> <tr> <td>74058 - 3-OHCOTININE</td> <td>74071 - NORNICOTINE</td> </tr> <tr> <td>74061 - ANABASINE</td> <td></td> </tr> </table>	74069 - NICOTINE	74063 - COTININE	74058 - 3-OHCOTININE	74071 - NORNICOTINE	74061 - ANABASINE	
74069 - NICOTINE	74063 - COTININE						
74058 - 3-OHCOTININE	74071 - NORNICOTINE						
74061 - ANABASINE							

91610 NMR, Lipoprofile
DT
92186 NMR Lipoprofile
NT

Specimen:															
Collect:	One LipoTube Also Acceptable Red Top														
Submit:	4 mL (Min:2 mL) Serum. Submit Refrigerated. Submit in a Standard Transport Tube.														
Special Handling:	Fast 12 hours prior to collection. Separate from cells within 8 hours.														
Rejection Criteria:	Lipemic Samples Non-fasting Specimen Plasma														
Stability:	Ambient: 1 Day(s); Refrigerated: 2 Week(s); Frozen: Unacceptable; Incubated: Unacceptable														
Methodology:	Quantitative Enzymatic; Quantitative Nuclear Magnetic Resonance Spectroscopy; Detergent Solubilization														
Performed:	Sun-Sat														
Reported:	5-9 Day(s)														
CPT Codes:	83704 80061														
Interpretive Data:	Please see report for interpretive data.														
Components:	<table border="0"> <tr> <td>92187 - LDL PARTICLE NUM</td> <td>92188 - SMALL LDL</td> </tr> <tr> <td>92189 - LARGE VLDL</td> <td>92190 - HDL PARTICLE NUM</td> </tr> <tr> <td>92191 - LARGE HDL</td> <td>92192 - LDL PARTICLE SIZE</td> </tr> <tr> <td>92193 - VLDL PARTICLE SIZE</td> <td>92194 - HDL PARTICLE SIZE</td> </tr> <tr> <td>92195 - TOTAL CHOLESTEROL</td> <td>92196 - TRIGLYCERIDES</td> </tr> <tr> <td>92197 - HDL CHOLESTEROL</td> <td>92198 - LDL CHOLESTEROL</td> </tr> <tr> <td>92199 - EER LIPOPROFILE</td> <td></td> </tr> </table>	92187 - LDL PARTICLE NUM	92188 - SMALL LDL	92189 - LARGE VLDL	92190 - HDL PARTICLE NUM	92191 - LARGE HDL	92192 - LDL PARTICLE SIZE	92193 - VLDL PARTICLE SIZE	92194 - HDL PARTICLE SIZE	92195 - TOTAL CHOLESTEROL	92196 - TRIGLYCERIDES	92197 - HDL CHOLESTEROL	92198 - LDL CHOLESTEROL	92199 - EER LIPOPROFILE	
92187 - LDL PARTICLE NUM	92188 - SMALL LDL														
92189 - LARGE VLDL	92190 - HDL PARTICLE NUM														
92191 - LARGE HDL	92192 - LDL PARTICLE SIZE														
92193 - VLDL PARTICLE SIZE	92194 - HDL PARTICLE SIZE														
92195 - TOTAL CHOLESTEROL	92196 - TRIGLYCERIDES														
92197 - HDL CHOLESTEROL	92198 - LDL CHOLESTEROL														
92199 - EER LIPOPROFILE															

70118 Opiates, Urine (Confirmation)
DT
74240 Opiates, Urine Confirmation
NT

Specimen:		
Collect:	Random Urine in Sterile Specimen Container	
Submit:	10 mL (Min:3 mL) Random Urine. Submit Refrigerated.	
Special Handling:	Avoid Repeated Freeze/Thaw Cycles	
Rejection Criteria:	Urine stabilized with additives or preservatives	
Stability:	Ambient: 1 Week(s); Refrigerated: 1 Week(s); Frozen: 24 Month(s); Incubated: Unacceptable	
Methodology:	Liquid Chromatography Mass Spectrometry	
Performed:	Mon-Fri	
Reported:	1-6 Day(s)	
CPT Codes:	80361	
Interpretive Data:	Please see report for interpretive data.	
Components:	74032 - CODEINE 74036 - HYDROCODONE 74049 - NORHYDROCODON 74121 - NOROXMORPHONE 74123 - TAPENTADOL-S 74126 - ZOLPIDEM	74043 - MORPHINE 74037 - HYDROMORPHONE 74119 - MEPERIDINE-U-Q 74122 - TAPENTADOL 74124 - TRAMADOL 74028 - 6ACETYLMORPH-U-Q

70420 Oxycodone, Urine (Confirmation)
DT
74250 Oxycodone, Urine Confirmation
NT

Specimen:		
Collect:	Random Urine in Sterile Specimen Container	
Submit:	10 mL (Min:3 mL) Random Urine. Submit Refrigerated.	
Special Handling:	Avoid Repeated Freeze/Thaw Cycles	
Rejection Criteria:	Urine stabilized with additives or preservatives	
Stability:	Ambient: 1 Week(s); Refrigerated: 1 Month(s); Frozen: 24 Month(s); Incubated: Unacceptable	
Methodology:	Liquid Chromatography Mass Spectrometry	
Performed:	Mon-Fri	
Reported:	1-6 Day(s)	
CPT Codes:	80365	
Interpretive Data:	Please see report for interpretive data.	
Components:	74052 - OXYCODONE 74051 - NOROXYCODONE	74053 - OXYMORPHONE

70126 Phencyclidine, Urine (Confirmation)
DT
74027 Phencyclidine, Urine Confirmation
NT

Specimen:	
Collect:	Random Urine in Sterile Specimen Container
Submit:	10 mL (Min:3 mL) Random Urine. Submit Refrigerated.
Special Handling:	Avoid Repeated Freeze/Thaw Cycles
Rejection Criteria:	Urine stabilized with additives or preservatives
Stability:	Ambient: 1 Week(s); Refrigerated: 1 Month(s); Frozen: 24 Month(s); Incubated: Unacceptable
Methodology:	Liquid Chromatography Mass Spectrometry
Performed:	Mon-Fri
Reported:	1-6 Day(s)
CPT Codes:	83992
Interpretive Data:	Please see report for interpretive data.

74054 Propoxyphene, Urine Confirmation
NT

Specimen:	
Collect:	Random Urine in Sterile Specimen Container
Submit:	10 mL (Min:3 mL) Random Urine. Submit Refrigerated.
Special Handling:	Avoid Repeated Freeze/Thaw Cycles
Rejection Criteria:	Urine stabilized with additives or preservatives
Stability:	Ambient: 1 Week(s); Refrigerated: 1 Month(s); Frozen: 24 Month(s); Incubated: Unacceptable
Methodology:	Liquid Chromatography Mass Spectrometry
Performed:	Mon-Fri
Reported:	1-6 Day(s)
CPT Codes:	80367
Interpretive Data:	Please see report for interpretive data.

