



## INTERPATH LABORATORY SERVICES

### TEST CHANGES EFFECTIVE NOVEMBER 15<sup>th</sup>, 2010

Please note the following test updates. Make the changes in your Interpath Service Manual if appropriate. Contact Client Services at 800-700-6891 for additional information.

<u>Test #</u>	<u>Test Name</u>	<u>Change/New Information</u>
2306	A-1-Antitrypsin	<b>New Reference Range</b> = 84-218 mg/dL
93628	AFFIRM	<b>Discontinued test.</b> See New Test # 1145
1145	AFFIRM	<b>New Test</b> Components: 1146 Candida 1147 Gardnerella 1148 Trichomonas Collect: Genital Swab Submit: In AFFIRM transport media ONLY Submit: Ambient Reference Range – Negative Stability: Ambient: 3 days Refrigerated: 3 days Frozen: Unacceptable CPT: 87480, 87510, 87660
2012	Amylase	<b>New Reference Range</b> = 29-103 U/L
2315	ASO Quant	<b>New Reference Range</b> = 0-250 IU/mL
91588	Calprotectin	<b>New Test</b> Collect Random Stool in clean container Submit 20 grams (Min 10 g) Ambient Reference Range: Normal = 0-50 ug/g Borderline = 51-120 ug/g Abnormal = 121 ug/g or greater Stability: Ambient: 5 days Refrigerated: 5 days Frozen: 1 year CPT: 83993

**UPDATE...UPDATE...UPDATE...UPDATE...UPDATE...UPDATE...UPDATE...UPDATE**

<b>2615</b>	<b>Ceruloplasmin</b>	<b>New Reference Range = 18-58 mg/dL</b>
<b>2356</b>	<b>Complement C3</b>	<b>New Reference Range</b> for Age 11 years or greater = 87-200 mg/dL
<b>2357</b>	<b>Complement C4</b>	<b>New Reference Range</b> for Age 11 years or greater = 19-52 mg/dL
<b>91359</b>	<b>Cortisol, Saliva</b>	<b>New Test</b> Collect Saliva between 2300-0100, using cotton Salivette. Rinse mouth thoroughly with water 10 minutes before collection. Place Salivette under tongue for 2 minutes, until swab is saturated. State collection time. Submit Frozen. Sample must be frozen within one hour of collection. Unacceptable if contaminated with blood, mucus, food or cellular debris. pH must be between 3.5-9.0. Reference Range: 0-0.111 ug/dL Stability: Ambient: Unacceptable Refrigerated: 8 hours Frozen: 6 months CPT: 82533
<b>2320</b>	<b>C-Reactive Protein</b>	<b>New Reference Range = 0-5 mg/L</b>
<b>2326</b>	<b>Immunoglobulin A</b>	<b>New Reference Range</b> for Ages 15 years or greater = 66-433 mg/dL
<b>2041</b>	<b>Lipase</b>	<b>New Reference Range = 11-82 U/L</b>
<b>2570</b>	<b>Lipoprotein (a)</b>	<b>Discontinued test.</b> See New Test # 9664
<b>9664</b>	<b>Lipoprotein (a)</b>	<b>New Test</b> Collect: One 7.5 mL SST tube Alternate: One 4mL Heparin or EDTA tube Submit: 1 mL (Min 0.5 mL) serum or plasma Submit: FROZEN Reference Range = 0-30 mg/dL Stability: Ambient: 8 hours Refrigerated:2 weeks Frozen: 3 months CPT: 83695

**UPDATE...UPDATE...UPDATE...UPDATE...UPDATE...UPDATE...UPDATE...UPDATE**

91572	Malaria Ab, IgG	<b>New Test</b> Collect: One 7.5 mL SST tube Submit: 1 mL serum (Min 0.5 mL) Submit: Refrigerated Avoid lipemia, hemolysis and icterus Reference Range = 0-1.00 IV Stability: Ambient: 2 days Refrigerated: 2 weeks Frozen: 1 year CPT: 86750
91167	Mumps, IgG	<b>Discontinued test.</b> See New Test # 2677
2677	Mumps, IgG	<b>New Test</b> Collect: One 7.5 mL SST tube Submit: 1 mL serum (Min 0.2 mL) Submit: Refrigerated Reference Range = >1.0 IV Stability: Ambient: 8 hours Refrigerated: 2 days Frozen: 1 month CPT: 86735
1012	Phosphorus, Inorganic	<b>New Reference Range =</b> Age: 0-10 years = 4.0-7.0 mg/dL Over 10 years = 2.5-5.0 mg/dL
94010	Rubeola, IgG	<b>Discontinued test.</b> See New Test # 2678
2678	Rubeola, IgG	<b>New Test</b> Collect: One 7.5 mL SST tube. Submit: 1 mL serum (Min 0.2 mL) Submit: Refrigerated Reference Range = >1.0 IV Stability: Ambient: 8 hours Refrigerated: 2 days Frozen: 1 month CPT:86765
1013	Uric Acid	<b>New Reference Range =</b> Male = 4.4-7.6 mg/dL Females = 2.3-6.6 mg/dL

**UPDATE...UPDATE...UPDATE...UPDATE...UPDATE...UPDATE...UPDATE...UPDATE**

<b>91228</b>	<b>Vitamin B6</b>	<b>New specimen requirement:</b> 1 ml (Min 0.5 mL) Heparin plasma Alternate: 1 mL serum (Do NOT use SST) CRITICAL Frozen and Protect from light. <b>New Reference Range</b> = 20-125 nmol/L
<b>91202</b>	<b>Very Long Chain Fatty Acid</b>	<b>Discontinued test.</b> Call Client Services at 800-700-6891 for new test information

**UPDATE...UPDATE...UPDATE...UPDATE...UPDATE...UPDATE...UPDATE...UPDATE**