



Interpath Laboratory, Inc. Test File Update

As an Interpath customer who receives electronic results or sends electronic orders you may need to be notified when we update our Service Manual. Although we try to keep these changes to a minimum, laboratory medicine is an evolving industry requiring changes to our technology from time to time. Depending on the requirements of your EMR or Hospital Information System you may be required to make similar changes to your system in order to correctly process inbound electronic results and create outbound electronic orders.

If you are uncertain that you are required to update your system we recommend that you contact your vendor for more information. As your laboratory service provider we are available to participate in the discussion with your vendor so that you clearly understand the impact of these changes.

Included in this email:

- This cover letter with a summary of changes
- Microsoft Word® Document with the detail of these changes to our Service Manual

Additional information including our most recent Service Manual and additional contact information can be found at www.interpathlab.com

Effective Date: May 20, 2021

Please note changes are to be implemented immediately.



Interpath Laboratory, Inc.

Test File Update

Order Code	Test Name	NC	CC	CPT	SRC	RRC	NT	DT	AOE
		Name Change	Component Change	CPT Change	Specimen Requirements Change	Reference Range Change	New Test	Discontinued Test	Ask on Order Entry Questions
74080	Buprenorphine, Urine Confirmation					◆			
74017	Cannabinoids, Urine Confirmation					◆			
74044	Naloxone, Urine Confirmation					◆			

74080 Buprenorphine, Urine Confirmation

Specimen:	
Collect:	Random Urine in Sterile Specimen Container
Submit:	10 mL (Min:3 mL) Random Urine. Submit Refrigerated.
Special Handling:	Avoid Repeated Freeze/Thaw Cycles
Rejection Criteria:	Urine stabilized with additives or preservatives
Stability:	Ambient: 1 Week(s); Refrigerated: 1 Month(s); Frozen: 24 Month(s); Incubated: Unacceptable
Methodology:	Liquid Chromatography Mass Spectrometry
Performed:	Mon-Fri
Reported:	1-6 Day(s)
CPT Codes:	80348
Interpretive Data:	Please see report for interpretive data.
Components:	74031 - BUPRENORPHINE 74047 - NORBUPRENORPH

Please note changes to buprenorphine cutoff.
Buprenorphine 5 ng/mL

74017 Cannabinoids, Urine Confirmation

Specimen:	
Collect:	Random Urine in Sterile Specimen Container
Submit:	10 mL (Min:3 mL) Random Urine. Submit Refrigerated.
Special Handling:	Avoid Repeated Freeze/Thaw Cycles
Rejection Criteria:	Urine stabilized with additives or preservatives
Stability:	Ambient: 1 Week(s); Refrigerated: 1 Month(s); Frozen: 24 Month(s); Incubated: Unacceptable
Methodology:	Liquid Chromatography Mass Spectrometry
Performed:	Mon-Fri
Reported:	1-6 Day(s)
CPT Codes:	80349
Interpretive Data:	Please see report for interpretive data.

Please note changes to THC Metabolite cutoff.
THC Metabolite (11N9COOHTHC) 15 ng/mL

74044 Naloxone, Urine Confirmation

Specimen:	
Collect:	Random Urine in Sterile Specimen Container
Submit:	10 mL (Min:3 mL) Random Urine. Submit Refrigerated.
Special Handling:	Avoid Repeated Freeze/Thaw Cycles
Rejection Criteria:	Urine stabilized with additives or preservatives
Stability:	Ambient: 1 Week(s); Refrigerated: 1 Month(s); Frozen: 24 Month(s); Incubated: Unacceptable
Methodology:	Liquid Chromatography Mass Spectrometry
Performed:	Mon-Fri
Reported:	1-6 Day(s)
CPT Codes:	80362
Interpretive Data:	Please see report for interpretive data.

Please note changes to Naloxone cutoff.
Naloxone 10 ng/mL