



INTERPATH
LABORATORY

Interpath Laboratory, Inc.
Test File Update

As an Interpath customer who receives electronic results or sends electronic orders you may need to be notified when we update our Service Manual. Although we try to keep these changes to a minimum, laboratory medicine is an evolving industry requiring changes to our technology from time to time. Depending on the requirements of your EMR or Hospital Information System you may be required to make similar changes to your system in order to correctly process inbound electronic results and create outbound electronic orders.

If you are uncertain that you are required to update your system we recommend that you contact your vendor for more information. As your laboratory service provider we are available to participate in the discussion with your vendor so that you clearly understand the impact of these changes.

Included in this email:

- This cover letter with a summary of the changes
- Microsoft Word® Document with the detail of these changes to our Service Manual
- Interpath Master Order/Result Compendium

Additional information including our most recent Service Manual and additional contact information can be found at www.interpathlab.com

Effective Date: December 15, 2025



Order Code	Test Name	NC	CC	CPT	SRC	RRC	NT	DT	AOE
		Name Change	Component Change	CPT Change	Specimen Requirements Change	Reference Range Change	New Test	Discontinued Test	Ask on Order Entry Questions
2073	Cortisol (AM)							◆	
2880	Cortisol						◆		
2081	Cortisol (PM)							◆	
2724	Cortisol, Random							◆	
91350	Culture, AFB [Blood, Bone Marrow]				◆				
2180	Dexamethasone Suppression Test							◆	
7020	Drug Screen, Urine without Confirmation		◆						
73030	Drug Screen, Targeted Urine Pain Management		◆						
93645	Heparin CoFactor II							◆	
74000	Medication Management Expanded Urine Confirmation		◆						
74005	Medication Management Expanded with TCA Urine Confirmation		◆						
2886	Post Dexamethasone						◆		
75001	THC-COOH, Serum					◆			



2073 Cortisol (AM)

DT

Please take note of Discontinued Test. If test is needed refer to new test 2880 Cortisol.

2880 Cortisol

NT

Specimen:	
Collect:	One SST Also Acceptable One Green Top (Li Heparin) One Lavender (EDTA) One Pink Top (EDTA) One Red Top One Standard Transport Tube
Submit:	1 mL (Min:0.5 mL) Serum in Standard Transport Tube. Submit Refrigerated. Submit in a Standard Transport Tube. Also Acceptable 1 mL (Min:0.5 mL) Serum in Standard Transport Tube. Submit Refrigerated. Submit in a Standard Transport Tube. 1 mL (Min:0.5 mL) Plasma in Standard Transport Tube. Submit Refrigerated. Submit in a Standard Transport Tube. 1 mL (Min:0.5 mL) Plasma in Standard Transport Tube. Submit Refrigerated. Submit in a Standard Transport Tube. 1 mL (Min:0.5 mL) Plasma in Standard Transport Tube. Submit Refrigerated. Submit in a Standard Transport Tube.
Special Handling:	Avoid Repeated Freeze/Thaw Cycles For patients receiving therapy with high biotin doses (>5 mg/day), no laboratory test specimen should be collected until at least 8 hours after the last biotin administration. Collection time: 6:00-10:00am. Reference Range: 6.02-18.40 ug/dL Collection time: 4:00-8:00pm. Reference Range: 2.68-10.50 ug/dL
Stability:	Ambient: 1 Day(s); Refrigerated: 4 Day(s); Frozen: 1 Year(s); Incubated: Unacceptable
Methodology:	Electrochemiluminescence Immunoassay (ECLIA)
Performed:	Mon-Fri
Reported:	1-3 Day(s)
CPT Codes:	82533

Please take note of New Test.

2081 Cortisol (PM)

DT

Please take note of Discontinued Test with no replacement.

2724 Cortisol, Random

DT

Please take note of Discontinued Test with no replacement.



91350 Culture, AFB [Blood, Bone Marrow]

NRC

Specimen:	
Collect:	Whole blood in Bactec Myco/F Lytic Also Acceptable Bone Marrow in Bactec Myco/F Lytic
Submit:	5 mL (Min:1 mL) Whole blood in Bactec Myco/F Lytic. Submit Ambient. Also Acceptable 5 mL (Min:1 mL) Bone Marrow in Bactec Myco/F Lytic. Submit Ambient.
Special Handling:	Critical Ambient
Stability:	Ambient: 3 Hour(s); Refrigerated: Unacceptable; Frozen: Unacceptable; Incubated: Unacceptable
Methodology:	BACTEC Continuous Monitoring System; Culture
Performed:	Sun-Sat
Reported:	Final: Negative at 43 days. Positive cultures are reported as soon as detected.
CPT Codes:	87116

Please take note of changes to Collect, Submit and Methodology.

2180 Dexamethasone Suppression Test

DT

Please take note of Discontinued Test with no replacement.

7020 Drug Screen, Urine without Confirmation

CC

Specimen:	
Collect:	Random Urine in Sterile Specimen Container Also Acceptable One Standard Transport Tube
Submit:	10 mL (Min:2 mL) Random Urine. Submit Refrigerated.
Special Handling:	Recommended pH is 3.0-11.0
Rejection Criteria:	Urine stabilized with additives or preservatives
Stability:	Ambient: 2 Hour(s); Refrigerated: 3 Day(s); Frozen: 12 Month(s); Incubated: Unacceptable
Methodology:	Homogenous Enzyme Immunoassay
Performed:	Sunday, Mon-Fri
Reported:	1-2 Day(s)
CPT Codes:	80307 82570
Interpretive Data:	Please see report for interpretive data.
Components:	7021 - AMPHETAMINES 7023 - BENZODIAZEPINES 7025 - COCAINE 7028 - OPIATES 7030 - OXYCODONE 7012 - CREATININE, URINE 7022 - BARBITURATES 7024 - CANNABINOIDS 7027 - METHADONE 7029 - PHENCYCLIDINE 7031 - ALCOHOL, URINE

Please take note of changes to Components.



73030 Drug Screen, Targeted Urine Pain Management

CC

Specimen:																																																						
Collect:	Random Urine in Sterile Specimen Container																																																					
	Also Acceptable One Standard Transport Tube																																																					
Submit:	10 mL (Min:4 mL) Random Urine. Submit Refrigerated.																																																					
Special Handling:	Avoid freeze and thaw cycles.																																																					
Stability:	Ambient: Unacceptable; Refrigerated: 3 Day(s); Frozen: 1 Week(s); Incubated: Unacceptable																																																					
Methodology:	Enzyme Immunoassay (EIA); Liquid Chromatography Mass Spectrometry																																																					
Performed:	Monday, Wednesday, Friday																																																					
Reported:	2-7 Day(s)																																																					
CPT Codes:	80307																																																					
Interpretive Data:	Please see report for interpretive data.																																																					
Components:	<table><tr><td>7022 - BARBITURATES</td><td>7031 - ALCOHOL, URINE</td></tr><tr><td>7024 - CANNABINOIDS</td><td>7012 - CREATININE, URINE</td></tr><tr><td>73021 - CARISOPRODOL</td><td>73011 - CLONAZEPAM</td></tr><tr><td>73039 - NORMEPERIDINE</td><td>73068 - METHYLPHENIDATE</td></tr><tr><td>73047 - NORBUPRENORPHINE</td><td>73048 - NORFENTANYL</td></tr><tr><td>73056 - MEPROBAMATE</td><td>73077 - DIHYDROCODEINE</td></tr><tr><td>73078 - PROPOXYPHENE</td><td>73012 - 6-ACETYLMORPHINE</td></tr><tr><td>73013 - 7-AMINOCLONAZEPAM</td><td>73014 - A-OH-ALPRAZOLAM</td></tr><tr><td>73016 - ALPRAZOLAM</td><td>73017 - AMPHETAMINE</td></tr><tr><td>73018 - A-OH-MIDAZOLAM</td><td>73019 - BENZOYLECGONINE</td></tr><tr><td>73022 - BUPRENORPHINE</td><td>73023 - CODEINE</td></tr><tr><td>73024 - DIAZEPAM</td><td>73026 - FENTANYL</td></tr><tr><td>73027 - GABAPENTIN</td><td>73028 - HYDROCODONE</td></tr><tr><td>73029 - HYDROMORPHONE</td><td>73032 - LORAZEPAM</td></tr><tr><td>73033 - MDA</td><td>73036 - MDEA-EVE</td></tr><tr><td>73037 - MDMA-ECSTASY</td><td>73042 - METHADONE</td></tr><tr><td>73043 - METHAMPHETAMINE</td><td>73044 - MIDAZOLAM</td></tr><tr><td>73049 - MORPHINE</td><td>73051 - NALOXONE</td></tr><tr><td>73052 - NORDIAZEPAM</td><td>73053 - NORHYDROCODONE</td></tr><tr><td>73058 - NOROXYCODONE</td><td>73059 - NOROXYMORPHONE</td></tr><tr><td>73061 - OXAZEPAM</td><td>73062 - OXYCODONE</td></tr><tr><td>73063 - OXYMORPHONE</td><td>73064 - PHENCYCLIDINE</td></tr><tr><td>73066 - PHENTERMINE</td><td>73067 - PREGABALIN</td></tr><tr><td>73069 - TAPENTADOL</td><td>73071 - TAPENTADOL-O-SULF</td></tr><tr><td>73072 - TEMAZEPAM</td><td>73073 - TRAMADOL</td></tr><tr><td>73074 - ZOLPIDEM</td><td>73076 - ZOLPIDEM P-4-COOH</td></tr></table>		7022 - BARBITURATES	7031 - ALCOHOL, URINE	7024 - CANNABINOIDS	7012 - CREATININE, URINE	73021 - CARISOPRODOL	73011 - CLONAZEPAM	73039 - NORMEPERIDINE	73068 - METHYLPHENIDATE	73047 - NORBUPRENORPHINE	73048 - NORFENTANYL	73056 - MEPROBAMATE	73077 - DIHYDROCODEINE	73078 - PROPOXYPHENE	73012 - 6-ACETYLMORPHINE	73013 - 7-AMINOCLONAZEPAM	73014 - A-OH-ALPRAZOLAM	73016 - ALPRAZOLAM	73017 - AMPHETAMINE	73018 - A-OH-MIDAZOLAM	73019 - BENZOYLECGONINE	73022 - BUPRENORPHINE	73023 - CODEINE	73024 - DIAZEPAM	73026 - FENTANYL	73027 - GABAPENTIN	73028 - HYDROCODONE	73029 - HYDROMORPHONE	73032 - LORAZEPAM	73033 - MDA	73036 - MDEA-EVE	73037 - MDMA-ECSTASY	73042 - METHADONE	73043 - METHAMPHETAMINE	73044 - MIDAZOLAM	73049 - MORPHINE	73051 - NALOXONE	73052 - NORDIAZEPAM	73053 - NORHYDROCODONE	73058 - NOROXYCODONE	73059 - NOROXYMORPHONE	73061 - OXAZEPAM	73062 - OXYCODONE	73063 - OXYMORPHONE	73064 - PHENCYCLIDINE	73066 - PHENTERMINE	73067 - PREGABALIN	73069 - TAPENTADOL	73071 - TAPENTADOL-O-SULF	73072 - TEMAZEPAM	73073 - TRAMADOL	73074 - ZOLPIDEM	73076 - ZOLPIDEM P-4-COOH
7022 - BARBITURATES	7031 - ALCOHOL, URINE																																																					
7024 - CANNABINOIDS	7012 - CREATININE, URINE																																																					
73021 - CARISOPRODOL	73011 - CLONAZEPAM																																																					
73039 - NORMEPERIDINE	73068 - METHYLPHENIDATE																																																					
73047 - NORBUPRENORPHINE	73048 - NORFENTANYL																																																					
73056 - MEPROBAMATE	73077 - DIHYDROCODEINE																																																					
73078 - PROPOXYPHENE	73012 - 6-ACETYLMORPHINE																																																					
73013 - 7-AMINOCLONAZEPAM	73014 - A-OH-ALPRAZOLAM																																																					
73016 - ALPRAZOLAM	73017 - AMPHETAMINE																																																					
73018 - A-OH-MIDAZOLAM	73019 - BENZOYLECGONINE																																																					
73022 - BUPRENORPHINE	73023 - CODEINE																																																					
73024 - DIAZEPAM	73026 - FENTANYL																																																					
73027 - GABAPENTIN	73028 - HYDROCODONE																																																					
73029 - HYDROMORPHONE	73032 - LORAZEPAM																																																					
73033 - MDA	73036 - MDEA-EVE																																																					
73037 - MDMA-ECSTASY	73042 - METHADONE																																																					
73043 - METHAMPHETAMINE	73044 - MIDAZOLAM																																																					
73049 - MORPHINE	73051 - NALOXONE																																																					
73052 - NORDIAZEPAM	73053 - NORHYDROCODONE																																																					
73058 - NOROXYCODONE	73059 - NOROXYMORPHONE																																																					
73061 - OXAZEPAM	73062 - OXYCODONE																																																					
73063 - OXYMORPHONE	73064 - PHENCYCLIDINE																																																					
73066 - PHENTERMINE	73067 - PREGABALIN																																																					
73069 - TAPENTADOL	73071 - TAPENTADOL-O-SULF																																																					
73072 - TEMAZEPAM	73073 - TRAMADOL																																																					
73074 - ZOLPIDEM	73076 - ZOLPIDEM P-4-COOH																																																					

Please take note of changes to Components.

93645 Heparin CoFactor II

DT

Please take note of Discontinued Test with NO replacement.



74000 Medication Management Expanded Urine Confirmation

CC

Specimen:		
Collect:	Random Urine in Sterile Specimen Container	
	Also Acceptable One Standard Transport Tube	
Submit:	10 mL (Min:3 mL) Random Urine. Submit Refrigerated.	
Special Handling:	Avoid Repeated Freeze/Thaw Cycles	
Rejection Criteria:	Urine stabilized with additives or preservatives	
Stability:	Ambient: Unacceptable; Refrigerated: 3 Day(s); Frozen: 1 Week(s); Incubated: Unacceptable	
Methodology:	Liquid Chromatography Mass Spectrometry	
Performed:	Mon-Fri	
Reported:	1-6 Day(s)	
CPT Codes:	82570	G0483
Interpretive Data:	Please see report for interpretive data.	
Components:	74066 - MDAMPHETAMINE 74124 - TRAMADOL 74054 - PROPOXYPHENE-U-Q 74051 - NOROXYCODONE 74053 - OXYMORPHONE 74036 - HYDROCODONE 74043 - MORPHINE 74119 - MEPERIDINE-U-Q 74048 - NORFENTANYL 74041 - METHADONE 74023 - MEPROBAMATE 74008 - 7AMCLONAZEPAM 74012 - LORAZEPAM 74014 - OXAZEPAM 74031 - BUPRENORPHINE 74004 - BUTALBITAL 74007 - SECOBARBITAL 74064 - METHAMPHET 74018 - ETHYLGLUCURON 74022 - CYCLOBENZAPRN-U-Q 74038 - LEVORPHANOL 74024 - KETAMINE 74042 - MITRAGYNE-U-Q 74072 - RITALINICACID 74029 - 6BNALTREXOL 74089 - NORTRIPTYLINE 74017 - 11N9COOHTHC-U-Q 74102 - ACNORFENTANYL 74067 - MDMETHAMPHET 74122 - TAPENTADOL 74027 - PHENCYCLIDINE-U-Q 74052 - OXYCODONE 74032 - CODEINE 74037 - HYDROMORPHONE 74049 - NORHYDROCODON 74034 - FENTANYL 74033 - METHADONEMET 74021 - CARISOPRODOL 74062 - COCAINEMET-U-Q 74009 - AOHALPRAZOLAM 74013 - NORDIAZEPAM 74016 - TEMAZEPAM 74047 - NORBUPRENORPH 74006 - PHENOBARBITAL 74059 - AMPHETAMINE 74028 - 6ACETYLMORPH-U-Q 74019 - ETHYLSULFATE 74003 - DEXMETHORPHAN 74001 - GABAPENTIN 74026 - NORKETAMINE 74068 - MPHENIDATE 74044 - NALOXONE-U-Q 74046 - NALTREXONE 74002 - PREGABALIN 74101 - ACFENTANYL 74126 - ZOLPIDEM-U-Q 7012 - CREATININE, URINE	

Please take note of changes to Components.



74005 Medication Management Expanded with TCA Urine Confirmation

CC

Specimen:

Collect:	Random Urine in Sterile Specimen Container																																																																	
	Also Acceptable One Standard Transport Tube																																																																	
Submit:	10 mL (Min:3 mL) Random Urine. Submit Refrigerated.																																																																	
Special Handling:	Avoid Repeated Freeze/Thaw Cycles																																																																	
Rejection Criteria:	Urine stabilized with additives or preservatives																																																																	
Stability:	Ambient: Unacceptable; Refrigerated: 3 Day(s); Frozen: 1 Week(s); Incubated: Unacceptable																																																																	
Methodology:	Liquid Chromatography Mass Spectrometry																																																																	
Performed:	Mon-Fri																																																																	
Reported:	1-6 Day(s)																																																																	
CPT Codes:	82570	G0483																																																																
Interpretive Data:	Please see report for interpretive data.																																																																	
Components:	<table><tr><td>74066 - MDAMPHETAMINE</td><td>74067 - MDMETHAMPHET</td></tr><tr><td>74124 - TRAMADOL</td><td>74122 - TAPENTADOL</td></tr><tr><td>74054 - PROPOXYPHENE-U-Q</td><td>74027 - PHENCYCLIDINE-U-Q</td></tr><tr><td>74051 - NOROXYCODONE</td><td>74052 - OXYCODONE</td></tr><tr><td>74053 - OXYMORPHONE</td><td>74032 - CODEINE</td></tr><tr><td>74036 - HYDROCODONE</td><td>74037 - HYDROMORPHONE</td></tr><tr><td>74043 - MORPHINE</td><td>74049 - NORHYDROCODON</td></tr><tr><td>74119 - MEPERIDINE-U-Q</td><td>74034 - FENTANYL</td></tr><tr><td>74048 - NORFENTANYL</td><td>74033 - METHADONEMET</td></tr><tr><td>74041 - METHADONE</td><td>74021 - CARISOPRODOL</td></tr><tr><td>74023 - MEPROBAMATE</td><td>74062 - COCAINEMET-U-Q</td></tr><tr><td>74008 - 7AMCLONAZEPAM</td><td>74009 - AOHALPRAZOLAM</td></tr><tr><td>74012 - LORAZEPAM</td><td>74013 - NORDIAZEPAM</td></tr><tr><td>74014 - OXAZEPAM</td><td>74016 - TEMAZEPAM</td></tr><tr><td>74031 - BUPRENORPHINE</td><td>74047 - NORBUPRENORPH</td></tr><tr><td>74004 - BUTALBITAL</td><td>74006 - PHENOBARBITAL</td></tr><tr><td>74007 - SECOBARBITAL</td><td>74059 - AMPHETAMINE</td></tr><tr><td>74064 - METHAMPHET</td><td>74028 - 6ACETYLMORPH-U-Q</td></tr><tr><td>74018 - ETHYLGLUCURON</td><td>74019 - ETHYLSULFATE</td></tr><tr><td>74022 - CYCLOBENZAPRN-U-Q</td><td>74003 - DEXMETHORPHAN</td></tr><tr><td>74038 - LEVORPHANOL</td><td>74001 - GABAPENTIN</td></tr><tr><td>74024 - KETAMINE</td><td>74026 - NORKETAMINE</td></tr><tr><td>74042 - MITRAGYNINE-U-Q</td><td>74068 - MPHENIDATE</td></tr><tr><td>74072 - RITALINICACID</td><td>74044 - NALOXONE-U-Q</td></tr><tr><td>74029 - 6BNALTREXOL</td><td>74046 - NALTREXONE</td></tr><tr><td>74089 - NORTRIPTYLINE</td><td>74002 - PREGABALIN</td></tr><tr><td>74073 - AMITRIPTYLINE</td><td>74077 - CLOMIPRAMINE</td></tr><tr><td>74078 - DESIPRAMINE</td><td>74081 - DOXEPIN</td></tr><tr><td>74086 - IMIPRAMINE</td><td>74079 - DESMETDOXEPIN</td></tr><tr><td>74094 - TRIMIPRAMINE</td><td>74017 - 11N9COOHTHC-U-Q</td></tr><tr><td>74101 - ACFENTANYL</td><td>74102 - ACNORFENTANYL</td></tr><tr><td>74126 - ZOLPIDEM-U-Q</td><td>7012 - CREATININE, URINE</td></tr></table>		74066 - MDAMPHETAMINE	74067 - MDMETHAMPHET	74124 - TRAMADOL	74122 - TAPENTADOL	74054 - PROPOXYPHENE-U-Q	74027 - PHENCYCLIDINE-U-Q	74051 - NOROXYCODONE	74052 - OXYCODONE	74053 - OXYMORPHONE	74032 - CODEINE	74036 - HYDROCODONE	74037 - HYDROMORPHONE	74043 - MORPHINE	74049 - NORHYDROCODON	74119 - MEPERIDINE-U-Q	74034 - FENTANYL	74048 - NORFENTANYL	74033 - METHADONEMET	74041 - METHADONE	74021 - CARISOPRODOL	74023 - MEPROBAMATE	74062 - COCAINEMET-U-Q	74008 - 7AMCLONAZEPAM	74009 - AOHALPRAZOLAM	74012 - LORAZEPAM	74013 - NORDIAZEPAM	74014 - OXAZEPAM	74016 - TEMAZEPAM	74031 - BUPRENORPHINE	74047 - NORBUPRENORPH	74004 - BUTALBITAL	74006 - PHENOBARBITAL	74007 - SECOBARBITAL	74059 - AMPHETAMINE	74064 - METHAMPHET	74028 - 6ACETYLMORPH-U-Q	74018 - ETHYLGLUCURON	74019 - ETHYLSULFATE	74022 - CYCLOBENZAPRN-U-Q	74003 - DEXMETHORPHAN	74038 - LEVORPHANOL	74001 - GABAPENTIN	74024 - KETAMINE	74026 - NORKETAMINE	74042 - MITRAGYNINE-U-Q	74068 - MPHENIDATE	74072 - RITALINICACID	74044 - NALOXONE-U-Q	74029 - 6BNALTREXOL	74046 - NALTREXONE	74089 - NORTRIPTYLINE	74002 - PREGABALIN	74073 - AMITRIPTYLINE	74077 - CLOMIPRAMINE	74078 - DESIPRAMINE	74081 - DOXEPIN	74086 - IMIPRAMINE	74079 - DESMETDOXEPIN	74094 - TRIMIPRAMINE	74017 - 11N9COOHTHC-U-Q	74101 - ACFENTANYL	74102 - ACNORFENTANYL	74126 - ZOLPIDEM-U-Q	7012 - CREATININE, URINE
74066 - MDAMPHETAMINE	74067 - MDMETHAMPHET																																																																	
74124 - TRAMADOL	74122 - TAPENTADOL																																																																	
74054 - PROPOXYPHENE-U-Q	74027 - PHENCYCLIDINE-U-Q																																																																	
74051 - NOROXYCODONE	74052 - OXYCODONE																																																																	
74053 - OXYMORPHONE	74032 - CODEINE																																																																	
74036 - HYDROCODONE	74037 - HYDROMORPHONE																																																																	
74043 - MORPHINE	74049 - NORHYDROCODON																																																																	
74119 - MEPERIDINE-U-Q	74034 - FENTANYL																																																																	
74048 - NORFENTANYL	74033 - METHADONEMET																																																																	
74041 - METHADONE	74021 - CARISOPRODOL																																																																	
74023 - MEPROBAMATE	74062 - COCAINEMET-U-Q																																																																	
74008 - 7AMCLONAZEPAM	74009 - AOHALPRAZOLAM																																																																	
74012 - LORAZEPAM	74013 - NORDIAZEPAM																																																																	
74014 - OXAZEPAM	74016 - TEMAZEPAM																																																																	
74031 - BUPRENORPHINE	74047 - NORBUPRENORPH																																																																	
74004 - BUTALBITAL	74006 - PHENOBARBITAL																																																																	
74007 - SECOBARBITAL	74059 - AMPHETAMINE																																																																	
74064 - METHAMPHET	74028 - 6ACETYLMORPH-U-Q																																																																	
74018 - ETHYLGLUCURON	74019 - ETHYLSULFATE																																																																	
74022 - CYCLOBENZAPRN-U-Q	74003 - DEXMETHORPHAN																																																																	
74038 - LEVORPHANOL	74001 - GABAPENTIN																																																																	
74024 - KETAMINE	74026 - NORKETAMINE																																																																	
74042 - MITRAGYNINE-U-Q	74068 - MPHENIDATE																																																																	
74072 - RITALINICACID	74044 - NALOXONE-U-Q																																																																	
74029 - 6BNALTREXOL	74046 - NALTREXONE																																																																	
74089 - NORTRIPTYLINE	74002 - PREGABALIN																																																																	
74073 - AMITRIPTYLINE	74077 - CLOMIPRAMINE																																																																	
74078 - DESIPRAMINE	74081 - DOXEPIN																																																																	
74086 - IMIPRAMINE	74079 - DESMETDOXEPIN																																																																	
74094 - TRIMIPRAMINE	74017 - 11N9COOHTHC-U-Q																																																																	
74101 - ACFENTANYL	74102 - ACNORFENTANYL																																																																	
74126 - ZOLPIDEM-U-Q	7012 - CREATININE, URINE																																																																	

Please take note of changes to Components.



2886 Post-Dexamethasone

NT

Specimen:	
Collect:	One SST Also Acceptable One Standard Transport Tube
Submit:	1 mL (Min:0.5 mL) Serum in Standard Transport Tube. Submit Refrigerated. Submit in a Standard Transport Tube.
Special Handling:	Avoid Repeated Freeze/Thaw Cycles For patients receiving therapy with high biotin doses (>5 mg/day), no laboratory test specimen should be collected until at least 8 hours after the last biotin administration. Minimize air exposure Reference range is applicable only to samples collected at 8:00am following administration of 1 mg dexamethasone at midnight.
Stability:	Ambient: 1 Day(s); Refrigerated: 4 Day(s); Frozen: 1 Year(s); Incubated: Unacceptable
Methodology:	Electrochemiluminescence Immunoassay (ECLIA)
Performed:	Mon-Fri
Reported:	1-3 Day(s)
CPT Codes:	82533

Please take note of New Test.

75001 THC-COOH, Serum

RRC

Specimen:	
Collect:	One Red Top Also Acceptable One Lavender (EDTA) One Pink Top (EDTA) One Standard Transport Tube
Submit:	2 mL (Min:1 mL) Serum in Standard Transport Tube. Submit Frozen. Submit in a Standard Transport Tube. Also Acceptable 2 mL (Min:1 mL) Plasma in Standard Transport Tube. Submit Frozen. Submit in a Standard Transport Tube. 2 mL (Min:1 mL) Plasma in Standard Transport Tube. Submit Frozen. Submit in a Standard Transport Tube.
Special Handling:	Separate from cells within 2 hours of collection
Rejection Criteria:	Hemolyzed specimens Use of separator tubes Whole blood Too low volume. Whole blood. Light blue (citrate) or yellow (SPS or ACD solution).
Stability:	Ambient: 1 Week(s); Refrigerated: 2 Week(s); Frozen: 12 Month(s); Incubated: Unacceptable
Methodology:	Quantitative Liquid Chromatography-Tandem Mass Spectrometry
Performed:	Mon-Fri
Reported:	2-8 Day(s)
CPT Codes:	80349
Interpretive Data:	General Reference Range : <10 ng/mL For medical purpose only. Methodology: Quantitative Liquid Chromatography Tandem Mass Spectrometry. Test developed and validated by Interpath Laboratory. Questions can be directed to laboratory.

Please take note of changes to Reference Range.