



INTERPATH
LABORATORY

Interpath Laboratory, Inc.
Test File Update

As an Interpath customer who receives electronic results or sends electronic orders you may need to be notified when we update our Service Manual. Although we try to keep these changes to a minimum, laboratory medicine is an evolving industry requiring changes to our technology from time to time. Depending on the requirements of your EMR or Hospital Information System you may be required to make similar changes to your system in order to correctly process inbound electronic results and create outbound electronic orders.

If you are uncertain that you are required to update your system we recommend that you contact your vendor for more information. As your laboratory service provider we are available to participate in the discussion with your vendor so that you clearly understand the impact of these changes.

Included in this email:

- This cover letter with a summary of the changes
- Microsoft Word® Document with the detail of these changes to our Service Manual
- Interpath Master Order/Result Compendium

Additional information including our most recent Service Manual and additional contact information can be found at www.interpathlab.com

Effective Date: October 1, 2026



INTERPATH
LABORATORY

Interpath Laboratory, Inc.
Test File Update

Order Code	Test Name	NC	CC	CPT	SRC	RRC	NT	DT	AOE
		Name Change	Component Change	CPT Change	Specimen Requirements Change	Reference Range Change	New Test	Discontinued Test	Ask on Order Entry Questions
3796	Semen Analysis, Complete						◆		
3641	Semen Analysis, Complete							◆	



3796 Semen Analysis, Complete

NT

Specimen:																	
Collect:	Semen in Sterile Specimen Container																
Submit:	Semen in Sterile Specimen Container. Submit Ambient.																
Special Handling:	Deliver to laboratory within 30 minutes, but do not reject until greater than 1 hour. Abstain from sexual intercourse or masturbation for 48 hours, but no longer than 7 days prior to collection. Collect the specimen by masturbation and ejaculate into a clean plastic container provided by the laboratory. If masturbation is not possible, please consult your doctor. Note: Reflexes to Vitality (test 3788) if Motility result is <40%. Additional charges will apply.																
Rejection Criteria:	Deliver sample within 30 minutes, will be rejected if greater than 1 hour Do not use a condom																
Stability:	Ambient: 1 Hour(s); Refrigerated: Unacceptable; Frozen: Unacceptable; Incubated: Unacceptable																
Methodology:	Computer Assisted Semen Analysis or Microscopy																
Performed:	Mon-Fri																
Reported:	1-2 Day(s) This test not performed at all locations contact Client Services at 800-700-6891 for possible locations.																
CPT Codes:	89320																
Interpretive Data:	Please see report for interpretive data.																
Components:	<table><tr><td>3665 - VISCOSITY</td><td>3667 - PH</td></tr><tr><td>3670 - VOLUME</td><td>3673 - ANALYSIS START</td></tr><tr><td>3676 - DAYS ABSTINENT</td><td>3679 - APPEARANCE</td></tr><tr><td>3681 - LIQUEFACTION</td><td>3689 - WBC CONC.</td></tr><tr><td>3691 - ROUND CELLS</td><td>3692 - SPERM COUNT</td></tr><tr><td>3693 - TOTAL MOTILITY</td><td>3694 - PROG. MOTILITY</td></tr><tr><td>3696 - NONPROG. MOTILITY</td><td>3697 - IMMOTILE</td></tr><tr><td>3698 - MORPHOLOGY-NORM</td><td></td></tr></table>	3665 - VISCOSITY	3667 - PH	3670 - VOLUME	3673 - ANALYSIS START	3676 - DAYS ABSTINENT	3679 - APPEARANCE	3681 - LIQUEFACTION	3689 - WBC CONC.	3691 - ROUND CELLS	3692 - SPERM COUNT	3693 - TOTAL MOTILITY	3694 - PROG. MOTILITY	3696 - NONPROG. MOTILITY	3697 - IMMOTILE	3698 - MORPHOLOGY-NORM	
3665 - VISCOSITY	3667 - PH																
3670 - VOLUME	3673 - ANALYSIS START																
3676 - DAYS ABSTINENT	3679 - APPEARANCE																
3681 - LIQUEFACTION	3689 - WBC CONC.																
3691 - ROUND CELLS	3692 - SPERM COUNT																
3693 - TOTAL MOTILITY	3694 - PROG. MOTILITY																
3696 - NONPROG. MOTILITY	3697 - IMMOTILE																
3698 - MORPHOLOGY-NORM																	

Please note this new test, this new test follows WHO 6th Edition guidelines.

3641 Semen Analysis, Complete

DT

Please note discontinued test. For alternative please see new test 3796 – Semen Analysis, Complete.